

KREATECH DIAGNOSTICS PRODUCT GUIDE





INTRODUCTION

DEAR VALUED CUSTOMER

Welcome to our new Product Guide 2010 ! This catalog includes a lot of new products which we have added into our portfolio throughout last year and we would like to specifically highlight the following sections:

- Our new line of CISH products (page 92)
- New multicolor probe panels for preimplantation genetic screening (page 66)
- Our unique target labeling products based on ULS[™]
- And many new FISH probes, e.g. for prostate cancer (TMPRSS2-ERG), glioma and sarcoma (CDK-4), our solid tumor probes MDM4, FGFR1, and our new MYC breakapart probe for lymphoma

We hope that you will find a lot of products in this catalog suitable for your daily research. If you have any additional questions do not hesitate to contact us or visit www.kreatech.com to find the most up-to-date information about our products.

Kreatech Diagnostics

PRODUCT SELECTION GUIDE



TABLE OF CONTENTS

Introduction	3
Table of Contents	5
Chromosome Index	6
HUGO Gene Symbols	8

MOLECULAR CYTOGENETICS

POSEIDON FISH Probes	10
Oncology Probes	10
Hematology	10
Solid Tumors	46
Preimplantation Genetic Screening	66
Prenatal Probes	68
Microdeletion Probes	70
Sub-Telomere Probes	82
Satellite Enumeration Probes	84
Whole Chromosome Probes	86
Arm Specific Probes	88
Mouse FISH Probes	90
Customized DNA Probes	91

CISH	92
Unistar	92
Twinstar	92

Cell Culture Media	94
Equipment	97
Kits and Reagents	98
FISH <i>Bright</i> Labeling Kits	100

MICROARRAY LABELING

The ULS Labeling Technology	103
DNA Labeling	104
RNA Labeling	106
Protein Labeling	108

General Nucleic Acid Labeling	109
Antibody Labeling	110
KREA <i>pure</i>	111

Technical Appendix	112
Terms and Conditions	114
CE marking & ISO certificates	115
Product Index	116
How to order	126
International Distributors	127



/	20pter	KBI-40236	83
	, SE20 D20Z1	KBI-20020 KBI-40062	84 66
	20q11	KBI-10203 KBI-10721 KBI-10733	27 61 48
	- PTPRT 20q12	KBI-10203	27
	- ZNF217 20q13	KBI-10733	48
	AURKA 20q13	KBI-10721	61
20 \	20qter	KBI-40237	83



MYC 8q24

8qter

8

44 50

KBI-40216 83

	/ SE15 D1524	KBI-20015 KBI-40063	84 66
	SNRPN 15q11 …	KBI-40109 KBI-40116	75 81
	🗅 UBE3A 15q11 🖤	KBI-40110 KBI-40116	76 81
	~ 15q22	KBI-10504 KBI-10508	38 38
	- PML 15q24	KBI-10302 KBI-40109 KBI-40110	30 75 76
	IGF1R 15q26	KBI-40116	81
15	15qter	KBI-40227	83



KBI-40238 83

21qter ····

SE13/21 KBI-20027 FOX01 13q14 KBI-10716	94 57
13q14 KBI-40001 KBI-40003 KBI-40005 KBI-40050 KBI-40050 KBI-40054	68 69 67 67 66
DLEU1 13q14	18 22

13qter

13

KBI-10502

KBI-40225 83

36

	— 7pter	KBI-40213	83
	– EGFR 7p11	KBI-10702	59
	- SE7 D7Z1	KBI-20007 KBI-20031 KBI-40061	84 84 66
	ELN 7q11	KBI-40111	76
	— 7q22	KBI-10202 KBI-10207 KBI-40111	26 27 76
	[_] MET 7q31	KBI-10719	60
	- 7q35	KBI-10202 KBI-10207	26 27
7	– 7qter	KBI-40214	83



1pter

CHD5 1p36

Cat.#

KBI-40201 83

KBI-10507 KBI-10712

Page

39 53



	9pter CDKN2A 9p21	KBI-40217 KBI-10402 KBI-10710	83 33 50
	SE9 9qh	KBI-20009 KBI-40062	84 66
	9q21	KBI-10402 KBI-10710	33 50
Ł	ABL1 9q34	KBI-10005 KBI-10006 KBI-10008 KBI-10009 KBI-10508	13 13 14 14 38
9	9qter	KBI-40218	83

PPARG 3p25 KBI-10707	51
SE3 D3Z1	84 66 21
EVI inv(3) 3q26 ··· KBI-10204	28
KBI-10205	28
TERC 3q26 KBI-10110	21
KBI-10704	50
BCL6 3q27 KBI-10607	44
KBI-10730	64
3qter KBI-40206	83

— 3pter

KBI-40205 83

CHROMOSOME INDEX



	12pter	KBI-40223	83
	ETV6 12p13	KBI-10401 KBI-10403	32 33
1	SE12 D12Z1	KBI-20012 KBI-40062	84 66
	GLI1 12q13	KBI-10104 KBI-10108	21 22
	CHOP 12q13	KBI-10714	56
	CDK4 12q13	KBI-10725	58
	MDM2 12q15	KBI-10717	58
12	12qter	KBI-40224	83



_	XYpter KBI-40240	83
	SEY DYZ3 KBI-20024 KBI-20030 KBI-20032	84 84 84
	KBI-40005 KBI-40008 KBI-40063	69 69 66
	Yq11 KBI-40051	67
-	SEY DYZ1 KBI-20025	84
Y	XYgter KBI-40241	83

CTNND2 5p15 ····· KBI-40106	78
SE1/5/19 KBI-20026	84 > 4
CDC25C 5q31 ····· KBI-10209	25
KDI-10210	20
EGRI 5q31 KBI-10208 KBI-10209	25
KBI-10210	26
KBI-10709	59
CSF1R 5q33 KBI-10209	25
KBI-10210	20
PDGFRB 5q33 KBI-10004 KBI-40064	17
RE14 Ec22 KBI 10200	25
KPS14 5q55 KBI-10209	25 26
NSD1 5a35	74
5ater KBI-40210	83
11pter KBI-40221	83
SE11 D11Z1 KBI-20011	84 66
NYEOV 11-12 KDI 1000E	25
MYEOV 11q13 KBI-10605	35
CCND1 11q13 KBI-10604 KBI-10609	$41 \\ 43$
KBI-10734	62
ATM 11q22 KBI-10103	19
KBI-10108	22
	25
	50
11q23 KBI-10502	20
MLL 11q23 KBI-10502 KBI-10502 KBI-10303 KBI-10711	30 53
MLL 11q23 KBI-10502 KBI-10502 KBI-10502 KBI-10502 KBI-10502 KBI-10502	30 53 83
MLL 11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222	30 53 83
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230	30 53 83 83
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40101	30 53 83 83 77 1
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40101 TP52 17p13 KBI-40101	30 53 83 77 15
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40101 KBI-10711 KBI-10011	30 53 83 77 15 19
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40101 TP53 17p13 KBI-10011 KBI-10113 KBI-10113	30 53 83 83 77 15 19 22
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40101 TP53 17p13 KBI-10111 KBI-10113 KBI-10114 KBI-10114 KBI-10114	30 53 83 83 77 15 19 22 23 37
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40101 TP53 17p13 KBI-10111 KBI-10113 KBI-10113 KBI-10178 KBI-10134	30 53 83 83 77 15 19 22 23 37 65
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40101 TP53 17p13 KBI-10111 KBI-10113 KBI-10113 KBI-1012 KBI-1013 KBI-10738 KBI-10738	30 53 83 83 77 15 19 22 23 37 65 61
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40101 TP53 17p13 KBI-10111 KBI-10112 KBI-10112 KBI-10113 KBI-10113 KBI-10738 KBI-10738 AURKB 17p13 KBI-10722 RAI1 17p11 KBI-40101	30 53 83 83 77 15 19 22 23 37 65 61 77
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40101 TP53 17p13 KBI-10011 KBI-10113 KBI-10114 KBI-10508 KBI-10738 AURKB 17p13 KBI-10722 RAI1 17p11 KBI-40101 SE17 D1721 KBI-20017	30 53 83 83 77 15 19 22 23 37 65 61 77 84
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40101 TP53 17p13 KBI-10111 KBI-10112 KBI-10112 KBI-10113 KBI-10114 KBI-10738 KBI-10722 AURKB 17p13 KBI-10722 RAI1 17p11 KBI-40101 SE17 D17Z1 KBI-20017 KBI-40062 KBI-40062	30 53 83 83 77 15 19 22 23 37 65 61 77 84 66
11q23 KBI-10302 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40101 TP53 17p13 KBI-10111 KBI-10113 KBI-10112 KBI-10114 KBI-10114 KBI-10738 AURKB 17p13 AURKB 17p13 KBI-10722 RAI1 17p11 KBI-40101 SE17 D17Z1 KBI-20017 KBI-4062 NF1 17q11	30 53 83 83 77 15 19 22 23 37 65 61 77 84 66 74
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40101 TP53 17p13 KBI-10111 KBI-10112 KBI-10112 KBI-10113 KBI-10114 KBI-1058 KBI-10738 AURKB 17p13 KBI-10722 RAI1 17p11 KBI-40101 SE17 D1721 KBI-20017 KBI-40062 NF1 17q11 NF1 17q12 KBI-40114 HER2 17q12 KBI-40114	30 53 83 83 77 15 19 22 23 37 65 61 77 84 66 74 47
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40101 TP53 17p13 KBI-10111 KBI-10112 KBI-10112 KBI-10113 KBI-10114 KBI-10738 KBI-10738 AURKB 17p13 KBI-10722 RAI1 17p11 KBI-40101 SE17 D17Z1 KBI-400172 NF1 17q11 KBI-40114 HER2 17q12 KBI-10701 KBI-10735 KBI-10735	30 53 83 83 77 15 19 22 23 37 65 61 77 84 66 74 47 48
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40101 TP53 17p13 KBI-10111 KBI-10112 KBI-10112 KBI-10113 KBI-10114 KBI-10738 KBI-10738 AURKB 17p13 KBI-10722 RAI1 17p11 KBI-40101 SE17 D17Z1 KBI-40017 KBI-40062 NF1 17q11 KBI-10735 KBI-10735 RARA 17q21 KBI-10302	30 53 83 83 77 15 19 22 23 37 65 65 61 77 84 66 74 47 48 30
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40101 TP53 17p13 KBI-10111 KBI-10112 KBI-10112 KBI-10113 KBI-10114 KBI-10738 KBI-10738 AURKB 17p13 KBI-10722 RAI1 17p11 KBI-40101 SE17 D17Z1 KBI-40017 KBI-40020 NF1 17q11 HER2 17q12 KBI-10701 KBI-10735 RARA 17q21 KBI-10305 KBI-10305	30 53 83 83 77 15 19 22 23 37 65 61 77 84 66 74 47 48 30 31
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40210 TP53 17p13 KBI-10111 KBI-10113 KBI-10112 KBI-10113 KBI-10113 KBI-10509 KBI-10738 AURKB 17p13 KBI-10722 RAI1 17p11 KBI-40101 SE17 D1721 KBI-40011 SE17 D1721 KBI-40021 KBI-10735 KBI-10735 RARA 17q21 KBI-10302 KBI-10325 KBI-10324 KBI-10724 KBI-10724	30 53 83 17 15 19 12 23 37 65 61 77 77 84 66 74 47 48 30 31 47 48
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40210 TP53 17p13 KBI-40101 TP53 17p13 KBI-10114 KBI-1012 KBI-10114 KBI-10735 KBI-10738 AURKB 17p13 KBI-10722 RAI1 17p11 KBI-40101 SE17 D17Z1 KBI-40021 KBI-40022 KBI-10735 RARA 17q21 KBI-10735 TOP2A 17q21 KBI-10735 MDD 17232 KDI 17224	30 53 83 17 15 19 122 23 37 65 61 77 77 18 466 16 74 17 47 18 30 11 47 18 47 18 47 18
11q23 KBI-10302 MLL 11q23 KBI-10303 KBI-10711 KBI-10303 11qter KBI-10303 KBI-10711 KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40101 TP53 17p13 KBI-10114 KBI-1012 KBI-10113 KBI-10735 KBI-10738 AURKB 17p13 KBI-10722 RAI1 17p11 KBI-40101 SE17 D1721 KBI-40062 NF1 17q11 KBI-10701 KBI-10735 KBI-10735 RARA 17q21 KBI-10305 TOP2A 17q22 KBI-10011 KBI-10735 KBI-10714 KBI-10735 KBI-10735	30 53 83 17 15 19 12 23 37 65 61 77 77 84 66 74 47 48 30 31 47 48 15 74
11q23 KBI-10302 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40101 TP53 17p13 KBI-10111 KBI-10112 KBI-10112 KBI-1012 KBI-10113 KBI-1013 KBI-10114 KBI-10735 KBI-10738 AURKB 17p13 KBI-10722 RAI1 17p11 KBI-40101 SE17 D17Z1 KBI-40021 KBI-10735 KBI-10735 RARA 17q21 KBI-10305 TOP2A 17q21 KBI-10735 MPO 17q22 KBI-10011 KBI-40101 KBI-4014	30 53 83 17 15 19 122 23 37 65 61 77 77 84 66 1 74 47 47 48 30 31 47 48 15 74 48 15 74 83
11q23 KBI-10502 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-40222 17pter KBI-40230 LIS 17p13 KBI-40210 TP53 17p13 KBI-10111 KBI-10112 KBI-10112 KBI-10113 KBI-10113 KBI-10114 KBI-10114 KBI-10735 KBI-10738 AURKB 17p13 KBI-10722 RAI1 17p11 KBI-40101 SE17 D1721 KBI-40021 KBI-10735 KBI-10735 RARA 17q21 KBI-10305 TOP2A 17q21 KBI-10735 MPO 17q22 KBI-10711 KBI-10735 KBI-10735	30 53 83 17 15 19 122 23 37 65 61 77 77 84 66 1 74 47 47 48 30 11 47 18 15 18 15 74 48 18 30 11 47 18 30 11 47 18 83 18
11q23 KBI-10302 MLL 11q23 KBI-10303 KBI-10711 KBI-10303 I1qter KBI-10711 11qter KBI-10222 17pter KBI-40230 LIS 17p13 KBI-10011 KBI-10113 KBI-10113 KBI-10113 KBI-10113 KBI-10114 KBI-10738 AURKB 17p13 KBI-10722 RAI1 17p11 KBI-40101 SE17 D1721 KBI-40101 SE17 D1721 KBI-40101 SE17 D1721 KBI-40104 HER2 17q12 KBI-10735 RARA 17q21 KBI-10724 KBI-10735 KBI-10724 MPO 17q22 KBI-10714 KBI-40114 17qter KBI-40231 Yynter	30 53 83 17 15 19 22 23 37 65 61 77 74 48 30 31 47 48 30 11 47 48 30 11 47 11 48 11 83 11
11q23 KBI-10302 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-10222 17pter KBI-40230 LIS 17p13 KBI-40101 TP53 17p13 KBI-10011 KBI-10113 KBI-10113 KBI-10113 KBI-10113 KBI-10113 KBI-10113 KBI-10113 KBI-10738 AURKB 17p13 KBI-10722 RAI1 17p11 KBI-40101 SE17 D1721 KBI-40101 SE17 D1721 KBI-107035 RARA 17q21 KBI-107035 MPO 17q22 KBI-10735 MPO 17q22 KBI-40114 17qter KBI-40240 XYpter KBI-40240	30 53 83 17 15 19 22 23 37 65 61 77 74 48 30 31 47 48 30 31 47 48 33 18 74 18 83 20
11q23 KBI-10302 MLL 11q23 KBI-10303 KBI-10711 KBI-10711 11qter KBI-10220 17pter KBI-40230 LIS 17p13 KBI-40101 TP53 17p13 KBI-10011 KBI-10112 KBI-10112 KBI-10113 KBI-10113 KBI-10113 KBI-10738 AURKB 17p13 KBI-10722 RAI1 17p11 KBI-40101 SE17 D1721 KBI-40101 SE17 D1721 KBI-40101 SE17 D1721 KBI-10701 KBI-10735 RARA 17q21 KBI-10735 KBI-10735 MPO 17q22 KBI-10714 KBI-40114 17qter KBI-40114 17qter KBI-40240 SH040240 SH040220 KBI-40240	30 53 83 12 15 19 122 23 37 65 61 77 74 48 30 31 47 48 30 31 47 48 33 15 74 18 30 31 47 48 30 31 47 83 83 79 20 20
11q23 KBI-10302 MLL 11q23 KBI-10303 KBI-10711 KBI-10303 I1qter KBI-10711 11qter KBI-40220 17pter KBI-40230 LIS 17p13 KBI-10011 TP53 17p13 KBI-10011 KBI-10113 KBI-10112 KBI-10113 KBI-10113 KBI-10113 KBI-10113 KBI-10113 KBI-10738 AURKB 17p13 KBI-10722 RAI1 17p11 KBI-40101 SE17 D1721 KBI-40101 SE17 D1721 KBI-10701 KBI-10702 KBI-10302 KBI-10735 RARA 17q21 KBI-10735 MPO 17q22 MPO 17q22 KBI-10011 KBI-40114 17qter KBI-40240 SHOX Xp22 SHOX Xp22 KBI-40112 STS Xp22 KBI-40115	30 53 83 12 15 19 122 23 377 15 15 61 77 84 66 24 74 48 30 31 47 48 15 74 83 15 74 83 83 79 80 24
11q23 KBI-10302 MLL 11q23 KBI-10303 KBI-10711 KBI-10303 I1qter KBI-10711 11qter KBI-40220 17pter KBI-40230 LIS 17p13 KBI-10011 KBI-10113 KBI-10112 KBI-10113 KBI-10113 KBI-10113 KBI-10113 KBI-10113 KBI-10113 KBI-10113 KBI-10113 KBI-10113 KBI-10738 AURKB 17p13 KBI-40101 SE17 D1721 KBI-40101 SE17 D1721 KBI-40101 SE17 D1721 KBI-10701 KBI-10702 KBI-10302 KBI-10735 RARA 17q21 KBI-10305 TOP2A 17q21 KBI-10735 MPO 17q22 KBI-10011 KBI-40114 17qter KBI-40240 SHOX Xp22 KBI-40115 SKAL Xp22 KBI-40115 KAI-40115	30 53 83 12 15 19 122 23 377 15 15 61 77 84 66 2 74 48 30 31 47 48 15 74 83 15 74 83 83 79 80 80

Page

83

Cat.#

KBI-40209

KBI-10208 KBI-10210 KBI-10709 KBI-40113

5pter

TERT 5p15

	PDGFRB 5q33	KBI-10004 KBI-40064	17 66
4	RPS14 5q33	KBI-10209 KBI-10210	25 26
_	NSD1 5q35	KBI-40113	74
E	5qter	KBI-40210	83
5			
	11pter	KBI-40221	83
	• SE11 D11Z1	KBI-20011 KBI-40062	84 66
	MYEOV 11q13	KBI-10605	35
	CCND1 11q13	KBI-10604 KBI-10609 KBI-10734	41 43 62
	ATM 11q22	KBI-10103 KBI-10108 KBI-10114	19 22 23
	11q23	KBI-10502	36
1	MLL 11q23	KBI-10303 KBI-10711	30 53
11	11qter	KBI-40222	83
/	17pter	KBI-40230	83
/	LIS 17p13	KBI-40101	77
	[,] TP53 17p13	KBI-10011 KBI-10112 KBI-10113 KBI-10114	15 19 22 23
		KBI-10509 KBI-10738	37 65
2	AURKE 17p13	KBI-10722	61 77
	SE17 D1771	KBI-90101	84
25		KBI-40062	66
	NF1 17q11	KBI-40114	74
	HER2 17q12	KBI-10701 KBI-10735	47 48
	RARA 17q21	KBI-10302 KBI-10305	30 31
17	TOP2A 17q21	KBI-10724 KBI-10735	47
//	MPO 17q22	KBI-10011 KBI-40114	15 74
	1/qter	кы-40231	83
	XYpter	KBI-40240	83
	SHUX XP22	KBI-40112	/9
	KAI 1 Xn22	KBI-40115	00 80
	TFE3 Xp11	KBI-10711	65
	•	KBI-40051	67
	SEX DXZ1	KBI-20023	84 84
		KBI-20032	84
		KBI-40003	69
	VICT	KBI-40063	66
	XIST Xq13	кы-40108	79
X	XYqter	KBI-40241	83

	4pter	KBI-40207	83
	FGFR3 4p16	KBI-10602	36
	WHSC1 4p16	KBI-40107	78
	SE4 D4Z1	KBI-20004 KBI-40061	84 66
	FIP1L1-CHIC2- PDGFRA 4q12	KBI-10003 KBI-10007	16 16
Ξ.			
	4qter	KBI-40208	83
4			

Cat.#

Page



16pter

FUS 16p11

SE16 D16Z2

CBFB 16q22 ...

MAF 16q23 -----

16qter

16

KBI-40228 83

KBI-10715 57

KBI-10304 31

KBI-10610 35

KBI-40229 83

84

67 66

KBI-20016

KBI-40050 KBI-40063



HUGO GENE SYMBOLS

The table below gives the approved HUGO Gene Symbol and aliases used for every gene in our Poseidon FISH DNA probes according to the HUGO Gene Nomenclature Committee (HGNC).

Chromosome position	Genes involved	HUGO Gene Symbols	Aliases	Cat.#	Page
1n36	CHD5	CHD5		KBI-10507	39
1032	MDM4	MDM4	ΜΡΜΥ	KBI-10712	53 54
1q21	S100A10	S100A10		KBI-10503	37
1q21 2n24	MYCN	MYCN		KBI-10507	39 52
2q11	LAF	AFF3	MLLT2-like, LAF4	KBI-10706	52
3p25	PPARy	PPARG	PPARG1, PPARG2, NR1C3, PPARgamma	KBI-10707	51
3q26	EVI	MECOM	MDS1-EVI1, PRDM3	KBI-10204	28
3q26	hTERC	TERC	TR, hTR, TRC3, SCARNA19	KBI-10110 KBI-10204 KBI-10205 KBI-10704	21 28 28 50
3q27	BCL6	BCL6	ZBTB27, LAZ3, BCL5, BCL6A	KBI-10607	44
4p16	FGFR3	FGFR3	CEK2, JTK4, CD333	KBI-10730	36
4p16	WHSC1	WHSC1	MMSET, NSD2	KBI-40107	78
4q12	CHIC2	CHIC2	BTL	KBI-10003	16
4q12	FIP1L1	FIP1L1	DKFZp586K0717	KBI-10003	16
4a12	DOCEDA	DDCERA		KBI-10007	16
4q12				KBI-10007	16
5015	CINND	CTNNDZ	NPKAP, 0124	KBI-40106 KBI-10208	78 25
5p15	hTERT	TERT	TRT. TP2. TCS1. hEST2. EST2	KBI-10210	26
-1			····, ·· =, · == ·, ··=· =, == · =	KBI-10709 KBI-40113	59 74
5.04	65 695 6	67 695 6		KBI-10209	25
5q31	CDC25C	CDC25C		KBI-10210 KBI-10709	26 59
				KBI-10208	25
5q31	EGR1		TIS8, G0S30, NGFI-A, KROX-24, ZIF-268, AT225, ZNF225	KBI-10209	25
				KBI-10709	59
5q33	CSF1R	CSF1R	C-FMS, CSFR, CD115	KBI-10209	25
F - 22	DDCEDD	DDCEDD		KBI-10210	17
5433	PDGFKB	PDGFKB	JTKTZ, CD1400, PDGFRT	KBI-40064	66
5q33	RPS14	RPS14	EMTB	KBI-10209 KBI-10210	25 26
5q35	NSD1	NSD1	ARA267, FLJ22263, KMT3B	KBI-40113	74
6a21	SEC63	SEC63	SEC63L. PRO2507. ERdi2	KBI-10105 KBI-10109	20 23
7.44		FOED		KBI-10504	38
/p11 7a11	EGFR, Her1 FLN	EGFR	ERBBT WBS_WS_SVAS	KBI-10702 KBI-40111	59 76
7q11	LIMK1	LIMK1	LIMK	KBI-40111	76
7q22	CUTL1	CUX1	CDP, CDP1, CUX, CUT, Clox, CDP/Cut, CDP/Cux, Cux/CDP, CASP, GOLIM6	KBI-10202	26 27
7q31	C-MET	MET	HGFR, RCCP2	KBI-10719	60
8p21	PNOC EGEP1	PNOC	PPNOC	KBI-10503	37
8q21	ETO	RUNX1T	CDR, ETO, MTG8, ZMYND2	KBI-10301	29
				KBI-10106	20
8q24	C-MYC	MYC	c-Myc, bHLHe39	KBI-10603	23 40
				KBI-10611	44
0.04	46	CDIVIDA		KBI-10704 KBI-10402	50 33
9p21	p16	CDKNZA	CDK41, p16, INK4a, M151, CMM2, ARF, p19, p14, INK4, p16INK4a, p19Arf	KBI-10710	50
				KBI-10005 KBI-10006	13
9q34	ABL	ABL1	JTK7, c-ABL, p150	KBI-10008	14
				KBI-10009 KBI-10508	14 38
				KBI-10508	38
9q34	ASS	ASS1	CTLN1	KBI-10006	13
				KBI-10009	14
10q23	PTEN	PTEN	MMAC1, TEP1, PTEN1	KBI-10718	60
11415	DCLI	CCNDT	021651	KBI-10604	41
11q13	CCND1	CCND1	U21B31	KBI-10605	35
				KBI-10609 KBI-10734	43 62
11q13	MYEOV	MYEOV	OCIM	KBI-10605	35
11a22	ATM	ATM	TEL1, TELO1	KBI-10103 KBI-10108	19 22
· ·			,	KBI-10114	23
11q23	MLL	MLL	TRX1, HRX, ALL-1, HTRX1, CXXC7, MLL1A, KMT2A	KBI-10303 KBI-10711	30 53

8

Chromosome position	Genes involved	HUGO Gene Symbols	Aliases	Cat.#	Page
11a23	PLZF	ZBTB16	PLZF	KBI-10502	36
12-12	TE1			KBI-10302	32
12µ13	IEL CDK4	EIVO		KBI-10403	33
12q13 12q13	CDK4 CHOP	CDK4 DDIT3	PSK-J3 CHOP10 GADD153 CHOP	KBI-10/25 KBI-1071/	58
12012	CII	CUI		KBI-10104	21
12413			UDM2 UDMY MCCE270	KBI-10108	22
12015	MDMZ	MDMZ	HDMZ, HDMX, MGC5370	KBI-10/17 KBI-10102	58 18
13q14	DLEU	DLEU1	LEU1, XTP6, NCRNA00021	KBI-10113	22
13a1/	FKHR	FOX01	EK H1	KBI-10502	36 57
15414	TKIIN	10/01	TKII	KBI-10601	24,39,42
				KBI-10602	36
				KBI-10603 KBI-10604	40 41
14q32	IGH	IGH@	IGHDY1, IGH	KBI-10605	35
				KBI-10606	41
				KBI-10729	30 63
15q11	SNRPN	SNRPN	SMN, SM-D, HCERN3, SNRNP-N, SNURF-SNRPN, RT-LI	KBI-40109	75
15q11	UBE3A	UBE3A	AS, ANCR, E6-AP, FLJ26981	KBI-40110 KBI-40116	76 81
15 a 2 2	SMADE	SMADE	U-T17400	KBI-10504	38
15422	SIVIADO	SIVIADO	RST17432	KBI-10508	38
15a24	PMI	PMI	MYL TRIM19 RNF71	KBI-10302 KBI-40109	30 75
15921				KBI-40110	76
15q26	IGF1R	IGF1R	JTK13, CD221, IGFIR, MGC18216, IGFR	KBI-40116	81
16g22	CBFB	CBFB	PEBP2B	KBI-10715	31
16q23	MAF	MAF	c-MAF	KBI-10610	35
17p13 17p13	AURKB	AURKB PAFAH1R1	Aik2, IPL1, AurB, AIM-1, ARK2, STK5 LIS1 ΡΔΕΔΗ	KBI-10722	61 77
17613	215			KBI-10011	15
				KBI-10112	19
17p13	p53	TP53	p53, LFS1	KBI-10113 KBI-10114	22
				KBI-10509	37
17p11	D A I 1	D A I 1	DKE7D4244120 SMCD SMS KIA41920 MGC12924	KBI-10738	65 77
17g11	NF1	NF1	DRI 2F434A133, SINCR, SINS, RIAA1020, MOC12024	KBI-40101	74
17a12	ERBB2	ERBB2	NEU, HER-2, CD340, HER2	KBI-10701	47
				KBI-10/35 KBI-10302	48 30
17q21	RARA	RARA	RAR, NR1B1	KBI-10305	31
17q21	TOP2A	TOP2A		KBI-10724	47
17~22	MDO	MDO		KBI-10733	15
17422	MPO	MPO	O/T	KBI-40114	74
18q11 18q21	SYI BCI2	SS18 BCL2	SYI Bcl-2	KBI-10/13 KBI-10606	56 41
18g21	MAIT	MALT1		KBI-10608	42
10921	MALI	MALTI		KBI-10731	64
19q13	CD37	CD37	TSPAN26	KBI-40064	66
20.44		1440054	594	KBI-10203	27
20011	MAPREI	MAPREI	EB1	KBI-10/21 KBI-10733	61 48
20q12	PTPRT	PTPRT	RPTPrho, KIAA0283	KBI-10203	27
20q13	AURKA	AURKA	BTAK, AurA, STK7, ARK1	KBI-10721	61
20q13 21g22	AML	RUNX1	PEBP2A2, AMLCR1	KBI-10733 KBI-10301	48 29
21q22	ERG	ERG	erg-3, p55	KBI-10726	63
21q22	TMPRSS2	TMPRSS2	PRSS10	KBI-10726	63 13
22~11	DCD	DCD	DOORED CMI DHI ALI	KBI-10005	13
22411	DCN	DCN	DZZ300Z, CML, FHL, ALL	KBI-10008	14
22a11	CI H22	CITCI 1		KBI-10009 KBI-40102	14 71 73
22a11	DGCR2	DGCR2	KIAA0163 LAN IDD DGS-C SE7-12	KBI-40102	71,73
22911	EWIS			KBI-40103	72,73
22412	L 4V 3			KBI-40102	71,73
22q13	SHANK3	SHANK3	SPANK-2, prosap2, KIAA1650, PSAP2	KBI-40103	72,73
Xn22	KAI	KAI1	KALIG-1	KBI-40104 KBI-40115	72,73 80
Xp22	SHOX	SHOX	PHOG, GCFX, SS, SHOXY	KBI-40112	79
Xp22	STS	STS	ARCS	KBI-40115	80 65
Xp11	TFE3	TFE3	TFEA, bHLHe33	KBI-10/41 KBI-40051	67
Xq13	XIST	XIST	NCRNA00001, DXS1089, swd66	KBI-40108	79

ONCOLOGY PROBES - HEMATOLOGY PROBES





From the 25,000 genes in the human genome, approximately 350 genes have been causally linked to the development of cancer. Variant or aberrant function of these so-called cancer genes may result from changes in genome copy number (through amplification, deletion, chromosome loss, or duplication), changes in gene and chromosome structure (through chromosomal translocation, inversion, or other rearrangements that lead to chimeric transcripts or deregulated gene expression) and point mutations (including base substitutions, deletions, or insertions in coding regions and splice sites). The vast majority (90%) of cancer genes are mutated or altered through chromosomal aberrations in somatic tissue, about 10% are altered in the germ line, thereby transmitting heritable cancer susceptibility through successive generations. In addition to high resolution chromosome banding and advanced chromosomal imaging technologies, chromosome aberrations in cancer cells can be analyzed with an increasing number of large-scale, comprehensive genomic and molecular genetic technologies – including fluorescence *in situ* hybridization (FISH).

Chromosomal translocation (t) is the process by which a break in at least two different chromosomes occurs, with exchange of genetic material between the chromosomes.

Description	Page
Chronic Myelogenous Leukemia (CML)	12
Chronic Lymphocytic Leukemia (CLL)	18
Myelodysplastic Syndrome (MDS)	24
Acute Myeloid Leukemia (AML)	29
Acute Lymphoblastic Leukemia (ALL)	32
Multiple Myeloma (MM)	34
Lymphoma	40

Translocation, Dual-Fusion Assay

Dual-fusion, dual-color FISH assays for translocation utilizes large probes that span 2 breakpoints or flanking regions on the different chromosomes. Dual-fusion, dual-color FISH is optimal for detection of low levels of nuclei possessing a simple balanced translocation, as it greatly reduces the number of normal background nuclei with an abnormal signal pattern.

Translocation, Break-Apart or Split Assay

FISH using dual-color, break-apart probes is very useful in the evaluation of genes known to have multiple translocation partners; the differently colored probes hybridize to targets on opposite sides of the breakpoint of the affected gene.

Translocation, Dual-Fusion Assay



Translocation, Break or Split Assay



Expected signal pattern:

In normal intact cells, two separate red and two separate green individual signals will be visible, whereas a reciprocal translocation will generate two fused red/green signals (often appearing as single yellow signals), accompanied by one red and one green signal (representing the normal chromosomes).

Expected signal pattern:

In normal cells two sets of red/green-fused signals (representing the two alleles) will be visible. In an abnormal diploid cell, in which one allele has been split by a translocation, a separated red and green signal will be visible in addition to the normal fused signal.

Chronic Myeloproliferative Disorders (CMPD)

Chromosomal translocations in chronic myeloproliferative diseases (CMPD) almost invariably result in expression of constitutively activated fusion tyrosine kinases. The hallmark of these diseases is CML, where the BCR/ABL activated tyrosine kinase results from the balanced reciprocal Philadelphia chromosome translocation t(9;22).

Chronic Myelogenous Leukemia (CML) - BCR/ABL t(9;22)

CML is a malignant chronic myeloproliferative disorder (MPD) of the hematopoietic stem cell. All CML have a t(9;22) causing fusion of the 3' ABL region at 9g34 with the 5' BCR region at 22q11. This chimeric BCR/ABL gene encodes a constitutively activated protein tyrosine kinase with profound effects on cell cycle, adhesion, and apoptosis. Understanding this process has led to the development of the drug imatinib mesylate (Gleevec[™]), the first in a new class of genetically targeted agents, this is a major advance in cancer treatment. Several different approaches are used to analyze the BCR/ABL t(9;22)(q34;q11) by FISH each providing different details about this translocation.

BCR/ABL Product Family

The Philadelphia chromosome is an abnormally short chromosome 22 that is one of the two chromosomes involved in a translocation with chromosome 9. This translocation t(9;22) (q34;q11) takes place in a single bone marrow cell and, through the process of clonal expansion, gives rise to the leukemia.

ABL and BCR are normal genes on chromosomes 9 and 22, respectively. The ABL gene encodes a tyrosine kinase enzyme whose activity is tightly controlled. In the formation of the Ph translocation, two fusion genes are generated: BCR-ABL on the Ph chromosome and ABL-BCR on the chromosome 9 participating in the translocation. The BCR-ABL gene encodes a protein with deregulated tyrosine kinase activity.

The presence of this protein in the CML cells is strong evidence of its pathogenetic role. The efficacy in CML of a drug that inhibits the BCRABL tyrosine kinase has provided the final proof that the BCR-ABL oncoprotein is the unique cause of CML. The Poseidon portfolio contains now 4 different probes for BCR/ABL to suit all needs for the detection of t(9;22) by FISH:

Interpretation guidelines for Poseidon **BCR/ABL** Probes





ON BCR/ABL t(9;22), Fusion

Cat.# KBI-10005 BCR/ABL t(9;22), Fusion

The BCR/ABL t(9;22) Fusion is optimized to detect the t(9;22) (q34;q11) reciprocal translocation in a dual-color, dual-fusion assay on metaphase/interphase spreads, blood smears and bone marrow cells.

This probe will also detect cryptic insertions of ABL into the BCR region not detectable by karyotyping and therefore described as Ph-negative.

ON BCR/ABL t(9;22), TC, D-Fusion

The BCR/ABL t(9;22), TC, D-Fusion probe is a triple-color, dualfusion probe build from the same regions as the dual-color, dual-fusion probe, but the proximal BCR region is labeled in blue. Using the triple-color probe allows to distinguish between the derivative chromosome 22, the Philadelphia chromosome, which will be observed as purple (red/blue) color, while the derivative chromosome 9 will show a yellow (red/green) signal.



BCR/ABL t(9;22) Fusion probe hybridized on patient material

showing t(9;22) (q34;q11) reciprocal translocation (2RG1R1G). Image kindly provided by Monika Conchon, São Paulo

Cat.# KBI-10006 BCR/ABL t(9;22), TC, D-Fusion





BCR/ABL t(9;22), TC, D-Fusion probe hybridized on patient material showing translocation of distal BCR (1BG1RB1R1G). Image kindly provided by Prof Siebert, Kiel

Morris et al, 1990, Blood, 76: 1812-1818. Dewald et al, 1998, Blood, 91: 3357-3365. Kolomietz et al, 2001, Blood, 97; 3581-3588. Huntly et al, 2003, Blood, 102; 1160-1168. Tkachuk et al., 1990, Science 250, 559-562

Ordering information	Color	Tests	Cat#
ON BCR/ABL t(9;22) Fusion	red/green	10	KBI-10005
ON BCR/ABL t(9;22) Fusion	red/green	20	KBI-12005

Literature:

Morris et al, 1990, Blood, 76: 1812-1818. Dewald et al, 1998, Blood, 91: 3357-3365. Kolomietz et al, 2001, Blood, 97; 3581-3588. Huntly et al, 2003, Blood, 102; 1160-1168. Tkachuk et al., 1990, Science 250, 559-562

Ordering information	Color	Tests	Cat#
ON BCR/ABL t(9;22) TC, D-Fusion	red/green/blue	10	KBI-10006

ON BCR/ABL t(9;22), DC, S-Fusion, ES

A single-fusion assay is preferably used for the initial screening of CML and ALL patients. By adding an additional region proximal to the breakpoints on chromosome 9q34, this probe will provide an extra signal on the der(9q34) in case of a t(9;22). The Philadelphia chromosome, der(22q), is visualized by the fusion signal.

ON BCR/ABL t(9;22), DC, S-Fusion

A simple dual-color, single-fusion assay is preferably used for the initial screening of CML and ALL patients. The Philadelphia chromosome, der(22q), is visualized by a fusion signal while the der(9q) shows no signal.

Cat.# KBI-10008 BCR/ABL t(9;22), Dual-Color, Single-Fusion, Extra Signal









BCR/ABL t(9;22), DC, S-Fusion, ES hybridized to a normal metaphase (2R2G).

Literature:

Morris et al, 1990, Blood, 76: 1812-1818. Dewald et al, 1998, Blood, 91: 3357-3365. Kolomietz et al, 2001, Blood, 97; 3581-3588. Huntly et al, 2003, Blood, 102; 1160-1168. Tkachuk et al., 1990, Science 250, 559-562

Ordering information	Color	Tests	Cat#
ON BCR/ABL t(9;22) DC, S-Fusion, ES	red/green	10	KBI-10008

Literature:

Morris et al, 1990, Blood, 76: 1812-1818. Dewald et al, 1998, Blood, 91: 3357-3365. Kolomietz et al, 2001, Blood, 97; 3581-3588. Huntly et al, 2003, Blood, 102; 1160-1168. Tkachuk et al., 1990, Science 250, 559-562

Ordering information	Color	Tests	Cat#
ON BCR/ABL t(9;22) DC, S-Fusion	red/green	10	KBI-10009

BCR/ABL t(9;22), DC, S-Fusion hybridized to a normal metaphase (2R2G).

CML secondary chromosomal changes

ON p53 (17p13) / MPO (17q22) "ISO 17q"

Isochromosome 17q is the most common isochromosome in cancer. It plays an important role in tumor development and progression. Hematologic malignancies such as chronic myeloid leukemia (CML) with isochromosome 17q carry a poor prognosis. Isochromosome 17q is the most common chromosome abnormality in primitive neuroectodermal tumors and medulloblastoma. Isochromosome 17q is, by convention, symbolized as i(17q).

The p53 (17p13) / MPO (17q22) "ISO 17q" probe is optimized to detect copy numbers of the p53 gene region at 17p13 and MPO gene region at 17q22. In case of i(17q) a signal pattern of three red signals for MPO (17q22) and one signal for p53 at 17p13 is expected.

Cat.# KBI-10011 p53 (17p13) / MPO (17q22) "ISO 17q"



p53 (17p13) / MPO (17q22) "ISO 17q" hybridized to a normal metaphase (2R2G).

Literature:

Becher et al, 1990, Blood, 75: 1679-1683. Fioretos et al, 1999, Blood, 94: 225-232.

Ordering information	Color	Tests	Cat#
ON p53 (17p13) / MPO (17q22) "ISO 17q"	green/red	10	KBI-10011

SE	8	(D8Z1)		
SE	7	(D7Z1)/	SE 8	(D8Z1)

Gain of chromosome 8 is the most common secondary chromosomal aberration in CML (approx. 34%).

Cat.# KBI-20008 SE 8 (D8Z1) Cat.# KBI-20031 SE 7 (D7Z1) / SE 8 (D8Z1)

See description under Satellite Enumeration probes on page 84.

Ordering information

SE 7 (D7Z1) / SE 8 (D8Z1)

SE 8 (D8Z1)

Color	Tests	Cat#
red/green	10	KBI-20008
red/green	10	KBI-20031

Other Myeloproliferative Diseases:

ON FIP1L1-CHIC2-PDGFRA (4q12) Del, Break

Idiopathic hypereosinophilic syndrome (HES) and chronic eosinophilia leukemia (CEL) represent the most recent additions to the list of molecularly defined chronic myeloproliferative disorders. A novel tyrosine kinase that is generated from fusion of the Fip1-like 1 (FIP1L1) and PDGFR α (PDGFRA) genes has been identified as a therapeutic target for imatinib mesylate in hypereosinophilic syndrome (HES). This fusion results from an approximately 800-kb interstitial chromosomal deletion that includes the cysteine-rich hydrophobic domain 2 (CHIC2) locus.

The FIP1L1-CHIC2-PDGFRA probe is optimized to detect the CHIC2 deletion at 4q12 associated with the FIP1L1/PDGFRA fusion in a Dual-Color, split assay. It also allows the detection of translocation involving the FIP1L1 and PDGFRA region. However, chromosome 4 polyploidy may provide additional signals not associated with a translocation involving 4q12.

ON FIP1L1-CHIC2-PDGFRA (4q12) Del, Break, TC

The FIP1L1-CHIC2-PDGFRA triple-color probe is optimized to detect the CHIC2 deletion at 4q12 associated with the FIP1L1/ PDGFRA fusion in a triple-color, split assay. It also allows the detection of translocation involving the FIP1L1 and PDGFRA region. The split of the green and blue signal will indicate a translocation at 4q12 independent of a possible chromosome 4 polyploidy.

Cat.# KBI-10003 FIP1L1-CHIC2-PDGFRA (4q12) Del, Break



FIP1L1-CHIC2-PDGFRA (4q12) Del, Break probe hybridized to a normal interphase/metaphase (2RG).

Literature:

Cools et al, N Engl J Med, 2003, 348, 1201-1214. Godlib et al, Blood, 2004, 103, 2879-2891.

Ordering information	Color	Tests	Cat#
ON FIP1L1-CHIC2-PDGFRA (4q12) Del, Break	red/green	10	KBI-10003

Cat.# KBI-10007 FIP1L1-CHIC2-PDGFRA (4q12) Del, Break, Triple-Color



FIP1L1-CHIC2-PDGFRA (4q12) Del, Break, TC probe hybridized to a normal metaphase (2BRG).

Literature:

Cools et al, 2003, N Engl J Med, 348: 1201-1214. Griffin et al, 2003, PNAS, 100: 7830-7835. Gotlib et al, 2004, Blood, 103; 2879-2891.

Ordering information	Color	Tests	Cat#
ON FIP1L1-CHIC2-PDGFRA (4q12) Del, Break, TC	red/green/blue	10	KBI-10007

ON PDGFRB (5q33) Break

PDGFRB activation has been observed in patients with chronic myelomonocytic leukemia/atypical chronic myeloid leukemia and has been associated with 11 translocation partners, the best known is the ETV6 or TEL gene on 12p13, causing a t(5;12) translocation. Cytogenetic responses are achieved with imatinib in patients with PDGFRB fusion positive, BCR/ABL negative CMPDs.

The PDGFRB probe is optimized to detect translocations involving the PDGFRB region at 5q33 in a dual-color, split assay.

ON FGFR1 (8p12) Break

FGFR1 has been implicated in the tumorigenesis of haematological malignancies, where it is frequently involved in balanced chromosomal translocations, including cases of chronic myeloid leukaemia (BCR-FGFR1 fusion) and the 8p11 myeloproliferative syndrome/stem cell leukaemia–lymphoma syndrome, which is characterized by myeloid hyperplasia and non-Hodgkin's lymphoma with chromosomal translocations fusing several genes, the most common being a fusion between ZNF198 and FGFR1.

Cat.# KBI-10004 PDGFRB (5q33) Break





Cat.# KBI-10737 FGFR1 (8p12) Break

See description under Solid Tumors on page 49.

PDGFRB (5q33) Break probe hybridized to a normal metaphase (2RG).

Literature:

Wlodarska et al, 1997, Blood, 89: 1716-1722. Wilkinson et al, 2003, Blood, 102: 4287-419

Ordering information	Color	Tests	Cat#
ON PDGFRB (5q33) Break	red/green	10	KBI-10004

Ordering information	Color	Tests	Cat#
FGFR1 (8p12) Break	red/green	10	KBI-10737

Chronic Lymphocytic Leukemia (CLL)

CLL accounts for about 30% of all leukemias in Europe and the USA. Distinct clonal chromosomal abnormalities can be identified in up to 90% of CLL cases of the B-cell lineage. By FISH the most common chromosomal changes in CLL and their frequencies have been identified as shown in the table below.

ON DLEU (13q14) / 13qter

Deletions of chromosome 13q14 have been reported not only in CLL but in a variety of human tumors, including other types of lymphoid and myeloid tumors, as well as prostate, head and neck, and non-small cell lung cancers. The deletion of 13q may be limited to a single locus (13q14), or accompanied with the loss of a larger interstitial region of the long arm of chromosome 13. A minimal critical region of 400 kb has been described containing the DLEU1, DLEU2 and RFP2 genes.

The DLEU (13q14) specific DNA probe is optimized to detect copy numbers of the DLEU gene region at 13q14. The 13qter (13q34) region is included to facilitate chromosome identification.

Del(13q14)	55%
Del(11q)	18%
Trisomy 12q	16%
Del(17p)	7%
Del(6q)	6%
Trisomy 8q	5%
t(14q32)	4%
Trisomy 3q	3%

Cat.# KBI-10102 DLEU (13q14) / 13qter



DLEU (13q14) / 13qter probe hybridized to patient material showing a 13q14 deletion (1R2G). Image kindly provided by Dr. Dastugue, Toulouse

Literature:

Wolf et al, 2001, Hum Mol Genet, 10: 1275-1285. Corcoran et al, 1998, Blood, 91: 1382-1390.

Ordering information	Color	Tests	Cat#
ON DLEU (13q14) / 13qter	red/green	10	KBI-10102

ON p53 (17p13) / SE 17

The p53 tumor suppressor gene at 17p13, has been shown to be implicated in the control of normal cellular proliferation, differentiation, and apoptosis. Allelic loss, usually by deletion, and inactivation of p53 have been reported in numerous tumor types and are associated with poor prognosis in CLL.

The p53 (17p13) specific DNA probe is optimized to detect copy numbers of the p53 gene region at 17p13. The chromosome 17 satellite enumeration probe (SE 17) at D17Z1 is included to facilitate chromosome identification.

ON ATM (11q22) / SE 11

Chromosome 11q22.3-23.1 deletions involving the ataxiateleangiectasia mutated (ATM) locus are detected at diagnosis in 15-20% of cases of B-cell chronic lymphocytic leukemia (CLL) and are associated with a relatively aggressive disease. Loss of the 11q22-23 region, where the ataxia-telangiectasia mutated (ATM) gene is located is also one of the most frequent secondary chromosomal aberrations in mantle cell lymphoma.

The ATM (11q22.3) specific DNA probe is optimized to detect copy numbers of the ATM gene region at region 11q22.3. The chromosome 11 satellite enumeration (SE 11) at D11Z1 probe is included to facilitate chromosome identification.

Cat.# KBI-10112 p53 (17p13) / SE 17





p53 (17p13) / SE 17 probe hybridized to a normal metaphase (2R2G).

Cat.# KBI-10103 ATM (11q22) / SE 11



ATM (11q22) / SE 11 hybridized to patient material showing a 11q22 deletion at the ATM gene region observed (1R2G). Image kindly provided by Dr Wenzel, Basel

Literature:

Amiel A et al, 1997, Cancer Gener.Cytogenet,, 97; 97-100 Drach J et al, 1998, Blood, 92; 802-809

Ordering information	Color	Tests	Cat#
ON p53 (17p13) / SE 17	red/green	10	KBI-10112
ON p53 (17p13) / SE 17	red/green	20	KBI-12112

Literature:

Doehner et al, 1997, Blood, 89: 2516-2522. Bigoni et al, 1997, Leukemia, 11: 1933-1940.

Ordering information	Color	Tests	Cat#
ON ATM (11q22) / SE 11	red/green	10	KBI-10103

ON 6q21 / SE 6

Deletions affecting the long arm of chromosome 6 (6q) are among the most commonly observed chromosomal aberrations in lymphoid malignancies and have been identified as adverse prognostic factor in subsets of tumors including CLL. A minimal deletion region has been identified within a 2-megabase segment of 6q21, between D6S447 and D6S246. The SEC63 gene is located within this critical region.

The 6q21 specific DNA probe is optimized to detect copy numbers of 6q at region 6q21.The chromosome 6 satellite enumeration probe (SE 6) at D6Z1 is included to facilitate chromosome identification.

ON C-MYC (8q24) / SE 8

The C-MYC gene produces an oncogenic transcription factor that affects diverse cellular processes involved in cell growth, cell proliferation, apoptosis and cellular metabolism. The C-MYC oncogene has been shown to be amplified in many types of human cancer such as bladder, breast and cervical. Amplification at 8q24 including C-MYC is also observed in 5% of CLL patients. C-MYC is also the prototype for oncogene activation by chromosomal translocation.

The C-MYC (8q24) specific DNA probe is optimized to detect copy numbers of the C-MYC gene region at 8q24. The chromosome 8 satellite enumeration probe (SE 8) at D8Z1 is included to facilitate chromosome identification.

Cat.# KBI-10105 6q21 / SE 6





6q21 / SE 6 probe hybridized to a normal metaphase (2R2G).

Literature:

Sherratt et al, 1997, Chromosome Res, 5: 118-124. Zhang et al, 2000, Genes Chrom Cancer, 27: 52-58.

Ordering information	Color	Tests	Cat#
ON 6q21 / SE 6	red/green	10	KBI-10105

Cat.# KBI-10106 C-MYC (8q24) / SE 8





C-MYC (8q24) / SE 8 hybridized to a normal metaphase (2R2G).

Literature:

Greil et al, 1991, Blood, 78: 180-191. Note: This probe should not be used to detect translocations involving C-MYC.

Ordering information	Color	Tests	Cat#
ON C-MYC (8q24) / SE 8	red/green	10	KBI-10106

ON hTERC (3q26) / 3q11

Amplification of the 3q26-q27 has a high prevalence in cervical, prostate, lung, and squamous cell carcinoma. This amplification can also be found to a lesser extent in CLL patients. The minimal region of amplification was refined within a 1- to 2-Mb genomic segment containing several potential cancer genes including hTERC, the human telomerase RNA gene.

The hTERC (3q26) specific DNA probe is optimized to detect copy numbers of the hTERC gene region at region 3q26. The 3q11 region probe is included to facilitate chromosome identification.

ON GLI (12q13) / SE 12

Trisomy 12 is the most common numerical chromosomal aberration in patients with B-cell chronic lymphocytic leukemia (B-CLL). Partial trisomy 12 of the long arm of chromosome 12 consistently includes a smaller region at 12q13-15 and has been observed in CLL and several other tumors. A number of loci located close to either MDM2 or CDK4/SAS, including the genes GADD153, **GLI**, RAP1B, A2MR, and IFNG, were found to be coamplified.

The GLI (12q13) specific DNA probe is optimized to detect copy numbers of the GLI gene region at region 12q13. The chromosome 12 satellite enumeration probe (SE 12) D12Z3 is included to facilitate chromosome identification.

Cat.# KBI-10110 hTERC (3q26) / 3q11





hTERC (3q26) / 3q11 probe hybridized to a normal interphase/ metaphase (2R2G).

Literature:

Arnold et al, 1996, Genes Chrom Cancer, 16: 46-54. Soder et al, 1997, Oncogene, 14: 1013-1021.

Ordering information	Color	Tests	Cat#
ON hTERC (3q26) / 3q11	red/green	10	KBI-10110

Cat.# KBI-10104 GLI (12q13) / SE 12





GLI (12q13) / SE 12 hybridized to patient material showing GLI (12q13) amplification (3R2G). Image kindly provided by Dr Wenzel, Basel

Literature:

Merup et al, 1997, Eur J Haematol, 58: 174-180. Dierlamm et al., 1997, Genes Chrom Cancer, 20: 155-166.

Ordering information	Color	Tests	Cat#
ON GLI (12q13) / SE 12	red/green	10	KBI-10104

CLL probe combinations

Most of the CLL probes are also available as combinations covering two critical loci in one hybridization. This is of particular interest if patient material is limited. The disadvantage is that no internal control is added. The following CLL probe combinations are available:

ON DLEU (13q14) / p53 (17p13)

Deletion of DLEU at 13q14 indicates a rather good prognosis, deletion of p53 at 17p13 is associated with poor prognosis.

ON ATM (11q22) / GLI (12q13)

Deletion of ATM at 11q22-q23 indicates a rather poor prognosis, amplification of GLI at 12q13 is associated with an intermediate prognosis.

Cat.# KBI-10113 DLEU (13q14) / p53 (17p13)





DLEU (13q14) / p53 (17p13) hybridized to a normal metaphase (2R2G).

Cat.# KBI-10108 ATM (11q22) / GLI (12q13)





ATM (11q22) / GLI (12q13) hybridized to a normal metaphase (2R2G).

Literature:

Amiel A et al, 1997, Cancer Gener.Cytogenet, 97; 97-100 Drach J et al, 1998, Blood, 92; 802-809 Stilgenbauer S et al, 1998, Oncogene, 16; 1891 – 1897 Wolf S et al, 2001, Hum. Molec. Genet., 10; 1275-1285

Ordering information	Color	Tests	Cat#
ON DLEU (13q14) / p53 (17p13)	red/green	10	KBI-10113

Literature:

Döhner H et al, 1997, Blood, 7; 2516-2522 Boultwood J, 2001, J. Clin. Pathol., 54; 512-516 Dierlamm J et al, 1998, Genes Chromosomes Cancer, 20; 155-166 Döhner H at al, 1999, J. Molec. Med., 77; 266-281

Ordering information	Color	Tests	Cat#
ON ATM (11q22) / GLI (12q13)	red/green	10	KBI-10108

ON 6q21 / MYC (8q24)

Deletion of 6q21 indicates an intermediate risk category, amplification of C-MYC at 8q24 is associated with poor prognosis.

ON p53 (17p13) / ATM (11q22)

Deletion of p53 and ATM are both indicating poor prognosis in CLL.

Cat.# KBI-10109 6q21 / MYC (8q24)



6q21 / MYC (8q24) hybridized to a normal interphase/ metaphase (2R2G).

Cat.# KBI-10114 p53 (17p13) / ATM (11q22)





p53 (17p13) / ATM (11q22) hybridized to a normal metaphase (2R2G).

Ordering information	Color	Tests	Cat#
ON 6q21 / C-MYC (8q24)	red/green	10	KBI-10109

Ordering information	Color	Tests	Cat#
ON p53 (17p13) / ATM (11q22)	red/green	10	KBI-10114

Other relevant CLL probes:

ON IGH (14q32) Break

Cat.# KBI-10601 IGH (14q32) Break

See description under Multiple Myeloma on page 39.

Ordering information	Color	Tests	Cat#
ON IGH (14q32) Break	red/green	10	KBI-10601

Myelodysplastic Syndromes (MDS)

The myelodysplastic syndromes (MDS) are a heterogeneous group of hematopoietic disorders characterized in most patients by peripheral blood cytopenia with hypercellular bone marrow and dysplasia of the cellular elements. Cytogenetic studies play a major role in confirmation of diagnosis and prediction of clinical outcome in MDS, and have contributed to the understanding of its pathogenesis. Clonal chromosomal abnormalities are detected by routine karyotyping techniques in 40%-70% of cases of de novo MDS, and 95% of cases of therapy-related MDS.

SE 12 (D12Z3)

Cat.# KBI-20012 SE 12 (D12Z3)

See description under Satellite Enumeration probes on page 85.

Ordering information	Color	Tests	Cat#
SE 12 (D12Z3)	red/green	10	KBI-20012

ON hTERT (5p15) / 5q31

The hTERT / 5q31 dual-color probe can be used to detect deletions involving band 5q31 in MDS and AML.

The 5q- specific DNA probe is optimized to detect copy numbers at the CDC25C/EGR1 gene region at 5q31. The hTERT gene region at 5p15 is included to facilitate chromosome identification.

ON MDS 5q- (5q31; 5q33)

The presence of del(5q), either as the sole karyotypic abnormality or as part of a more complex karyotype, has distinct clinical implications for myelodysplastic syndromes (MDS) and acute myeloid leukemia. Interstitial 5q deletions are the most frequent chromosomal abnormalities in MDS and are present in 10% to 15% of MDS patients. Two different critical regions are described, one at 5q31-q33 containing the CSF1R and RPS14 gene regions ,characteristic for the '5q-' syndrome, and a more proximal located region at 5q13-q31 containing the CDC25C and EGR1 gene regions.

The 5q- specific DNA probe is optimized to detect copy numbers at the CDC25C/EGR1 gene region at 5q31 and the CSF1R/RPS14 gene region at 5q33 simultaneously in a dual-color assay.

Cat.# KBI-10208 hTERT (5p15) / 5q31



hTERT (5p15) / 5q31 probe hybridized to a normal interphase/ metaphase (2R2G).

Literature:

Zhao et al, 1997, PNAS, 94; 6948-6053 Horrigan et al, 2000, Blood, 95; 2372-2377

Ordering information	Color	Tests	Cat#
ON hTERT (5p15) / 5q31	red/green	10	KBI-10208







MDS 5q- (5q31; 5q33) probe hybridized to patient material showing a 5q33 deletion (1R2G). Image kindly provided by Dr Mohr, Dresden

Literature:

Boultwood J e.a., Blood 2002; 99: 4638-4641 Zhao N e.a., PNAS 1997; 94: 6948-6953 Wang e.a., Haematologica 2008; 93: 994-1000 Ebert BL e.a., Nature 2008: 451: 335-339 Mohamedali A and Mufti GJ, Brit J Haematol 2008; 144: 157-168

Ordering information	Color	Tests	Cat#
ON MDS 5q- (5q31; 5q33)	red/green	10	KBI-10209

ON MDS 5q- (5q31; 5q33) / hTERT (5p15) TC

The 5q- specific DNA probe is optimized to detect copy numbers at the CDC25C/EGR1 gene region at 5q31 and the CSF1R/RPS14 gene region at 5q33 simultaneously in a dualcolor assay. The triple-color probe provides in addition to the two critical regions a control in blue targeting the hTERT gene region at 5p15.

ON MDS 7q- (7q22; 7q35)

Loss of a whole chromosome 7 or a deletion of the long arm, del(7q), are recurring abnormalities in malignant myeloid diseases. Most deletions are interstitial and there are two distinct deleted segments of 7q. The majority of patients have proximal breakpoints in bands q11-22 and distal breakpoints in q31-36. The CCAAT displacement protein (CUTL1) gene region is located in the 7q22 critical region.

The 7q- specific DNA probe is optimized to detect copy number of 7q at 7q22 and at 7q35 simultaneously in a dualcolor assay





MDS 5q- (5q31; 5q33) / hTERT (5p15) probe hybridized to a

normal metaphase (2R2G2B).

Color

red/green/blue

Tests Cat#

KBI-10210

10

Cat.# KBI-10202 MDS 7q- (7q22; 7q35)





MDS 7q- (7q22; 7q35) hybridized to patient material showing a 7q35 deletion (1R2G). Image kindly provided by Prof Jauch, Heidelberg

Literature:

Ordering information

ON MDS 5q- (5q31; 5q33) / hTERT (5p15) TC

Boultwood J e.a., Blood 2002; 99: 4638-4641 Zhao N e.a., PNAS 1997; 94: 6948-6953 Wang e.a., Haematologica 2008; 93: 994-1000 Ebert BL e.a., Nature 2008: 451: 335-339 Mohamedali A and Mufti GJ, Brit J Haematol 2008; 144: 157-168

Literature:

LeBeau et al., 1996, Blood, 88: 1930-1935. Doehner et al, 1998, Blood, 92: 4031-4035.

Ordering information	Color	Tests	Cat#
ON MDS 7q- (7q22; 7q35)	red/green	10	KBI-10202

5015

ON MDS 7q- (7q22; 7q35) / SE 7 TC

The 7q- specific DNA probe is optimized to detect copy number of 7q at 7q22 and at 7q35 simultaneously in a dual-color assay.

The chromosome 7 satellite enumeration probe (SE 7) at D7Z1 in blue is included to facilitate chromosome identification.

ON MDS 20q- (PTPRT 20q12) / 20q11

Acquired deletions of the long arm of chromosome 20 are found in several hematologic conditions and particularly in the myeloproliferative disorders (MPD) and myelodysplastic syndromes and acute myeloid leukemia (MDS/AML). A minimal critical region deleted in MPD and MDS has been identified at 20q12 which includes a protein tyrosine phosphatase receptor gene (PTPRT).

The 20q- (PTPRT, 20q12) specific DNA probe is optimized to detect copy numbers of 20q at region 20q12. A 20q11 region specific probe is included to facilitate chromosome identification.

Cat.# KBI-10207 MDS 7q (7q22; 7q35) / SE 7, Triple-Color





MDS 7q (7q22; 7q35) / SE 7 TC probe hybridized to a normal metaphase (2R2G2B).

Literature:

LeBeau et al., 1996, Blood, 88: 1930-1935. Doehner et al, 1998, Blood, 92: 4031-4035.







MDS 20q- (PTPRT 20q12) / 20q11 probe hybridized to a normal metaphase (2R2G).

Literature:

Bench et al, 2000, Oncogene, 19: 3902-3913. Asimakopoulos et al, 1994, Blood, 84: 3086-3094.

Tests	Cat#	Ordering information	Color
10	KBI-10207	ON 20q- (PTPRT 20q12) / 20q11	red/green

Tests Cat#

KBI-10203

10

ON EVI t(3;3); inv(3) (3q26) Break

The inv(3)(g21;g26) is a recurrent cytogenetic aberration of myeloid malignancy associated with fusion of EVI1 and RPN1 and a poor disease prognosis. Genomic breakpoints in 3q26 are usually located proximal to the EVI1 locus, spanning a region of several hundred kilobases. Other recurrent and sporadic rearrangements of 3q26 also cause transcriptional activation of EVI1 including the translocations t(3;3)(g21;g26) and t(3;21) (q26;q22). Breakpoints in the latter rearrangements span a wider genomic region of over 1 megabase encompassing sequences distal to EVI1 and neighboring gene MDS1.

The EVI t(3;3) inv(3) Break, dual-color probe is optimized to detect the inversion of chromosome 3 involving the EVI1 gene region at 3q26 in a dual-color, split assay on metaphase/ interphase spreads, blood smears and bone marrow cells.

Cat.# KBI-10204 EVI t(3;3); inv(3) (3q26) Break





EVI t(3;3);inv(3) (3q26) Break probe hybridized to patient material showing a rearrangement involving the EVI gene region at 3q26 (1RG1R1G). Image kindly provided by Dr. Reed, London

Literature[.]

Levy et al, 1994, Blood, 83: 1348-1354. Wieser et al, 2003, Haematologica, 88: 25-30. Melo et al, 2007, Leukemia, 22, 434-437.

Ordering information	Color	Tests	Cat#
ON EVI t(3;3); inv(3) (3q26) Break	red/green	10	KBI-10204

ON EVI t(3;3); inv(3), EVI (3q26) Break, TC

The EVI t(3;3) inv(3) Break, triple-color probe is optimized to detect the inversion of chromosome 3 involving the EVI gene region at 3q26 in a dual-color, split assay on metaphase/ interphase spreads, blood smears and bone marrow cells. By using a third color breakpoint variations can now be easily discovered.

Cat.# KBI-10205 EVI t(3;3); inv(3) (3q26), Triple-Color





EVI t(3;3); inv(3) (3q26) TC probe hybridized to patient material showing a rearrangement involving the EVI gene region at 3q26 (1RGB1B1RG). Image kindly provided by Prof Jauch, Heidelberg

Literature:

Levy et al, 1994, Blood, 83: 1348-1354. Wieser et al, 2003, Haematologica, 88: 25-30. Melo et al, 2007, Leukemia, 22, 434-437.

Ordering information	Color	Tests	Cat#
ON EVI t(3;3); inv(3) (3q26) Break, TC	red/green/blue	10	KBI-10205

Oncology Probes - Hematology Probes

Acute Myeloid Leukemia (AML)

Atleast 80% of patients with acute myeloid leukemia (AML) have an abnormal karyotype. Cytogenetic analysis provides some of the strongest prognostic information available, predicting outcome of both remission induction and postremission therapy. Abnormalities which indicate a good prognosis include t(8;21), inv(16), and t(15;17). Patients with AML that is characterized by deletions of the long arms or monosomies of chromosomes 5 or 7; by translocations or inversions of chromosome 3, t(6;9), t(9;22); or by abnormalities of chromosome 11q23 have particularly poor prognoses with chemotherapy.

ON AML/ETO t(8;21) Fusion

t(8;21)(q22;q22) is the most frequently observed karyotypic abnormality associated with acute myeloid leukemia (AML), especially in FAB M2. As a consequence of the translocation the AML1 (CBFA2, RUNX1) gene in the 21q22 region is fused to the ETO (MTG8, RUNX1T) gene in the 8q22 region, resulting in one transcriptionally active gene on the 8qderivative chromosome.

The AML/ETO t(8;21)(q21;q22) specific DNA probe is optimized to detect the reciprocal translocation t(8;21) in a dual-color, dual-fusion assay.

Cat.# KBI-10301 AML/ETO t(8;21) Fusion





AMI/ETO t(8;21) Fusion probe hybridized to a normal metaphase (2R2G).

Literature:

Sacchi et al, 1995, Genes Chrom Cancer, 79: 97-103. Hagemeijer et al, 1998, Leukemia, 12: 96-101.

Ordering information	Color	Tests	Cat#
ON AML/ETO t(8;21) Fusion	red/green	10	KBI-10301

ON PML/RARA t(15;17) Fusion

A structural rearrangement involving chromosomes 15 and 17 in acute promyelocytic leukemia (APL) was first recognized in 1977. The critical junction is located on the der(15) chromosome and consists of the 5' portion of PML fused to virtually all of the RARA gene. The PML/RARA fusion protein interacts with a complex of molecules known as nuclear co-repressors and histone deacetylase. This complex binds to the fusion protein and blocks the transcription of target genes. Other less common variant translocations fuse the RARA gene on 17q21 to the PLZF, NPM, NUMA, and STAT5b genes, respectively.

The PML/RARA t(15;17) specific DNA probe is optimized to detect the reciprocal translocation t(15;17) (q24;q21) in a dual-color, dual-fusion assay.

ON MLL (11q23) Break

The human chromosome band 11q23 is associated with a high number of recurrent chromosomal abnormalities including translocations, insertions, and deletions. It is involved in over 20% of acute leukemias. The MLL (Myeloid-Lymphoid Leukemia or Mixed-Lineage Leukemia) gene, named for its involvement in myeloid (usually monoblastic) and lymphoblastic leukemia, and less commonly in lymphoma is located in the 11q23 breakpoint region. Leukemias involving the MLL gene usually have a poor prognosis.

The MLL (11q23) break probe is optimized to detect translocations involving the MLL gene region at 11q23 in a dual-color split assay.



PML/RARA t(15,17) Fusion probe hybridized to a normal metaphase (2R2G).

Color

red/green

red/areen

Tests Cat#

10

20

KBI-10302

KBI-12302

Literature:

Ordering information

ON PML/RARA t(15,17) Fusion

ON PML/RARA t(15.17) Fusion

Schad et al, 1994, Mayo Clin Proc, 69: 1047-1053. Brockman et al, 2003, Cancer Genet Cytogenet, 145: 144-151.

Literature:

Kobayashi et al, 1993, Blood, 81: 3027-3022 Martinez-Climent et al, 1995, Leukemia, 9: 1299-1304.

Ordering information	Color	Tests	Cat#
ON MLL (11q23) Break	red/green	10	KBI-10303

MLL (11q23) Break probe hybridized to patient material

showing a translocation at 11q23 (1RG1R1G).

ON RARA (17q21) Break

This break apart probe can detect the numerous types of recurrent rearrangement of the RAR α (Retinoid acid receptor, alpha) gene with various gene partners (e.g., PML, NPM, MLL, FIP1L1, NuMA1, PLZF, amongst the others), leading to the formation of different reciprocal fusion proteins. The importance of retinoid metabolism in acute promyelocytic leukemia (APL) is highlighted by the numerous recent studies, but the different leukemogenic functions of the RAR α fusion proteins in the neoplastic myeloid development still has to be defined, as well as the distinct clinical outcome of the patients with the variant forms of APL.

ON CBFB t(16;16); inv(16) Break

Cat.# KBI-10304 CBFB t(16;16); inv(16) Break

Inv(16)(p13;q22) and t(16;16)(p13;q22) are recurring chromosomal rearrangements in AML. In both the inversion and translocation, the critical genetic event is the fusion of the CBFB gene at 16q22 to the smooth muscle myosin heavy chain (MYH11) at 16p13. A deletion of between 150 and 350 kb centromeric to the p-arm inversion breakpoint cluster region can be observed in some patients containing the 5' portion of the myosin heavy chain (MYH11) gene.

The CBFB t(16;16) inv(16) break probe is optimized to detect the inversion of chromosome 16 involving the CBFB gene region at 16q22 in a dual-color, split assay.



Ordering information	Color	Tests	Cat#
ON RARA (17q21) Break	red/green	10	KBI-10305

Ordering information	Color	Tests	Cat#
ON CBFB t(16;16); inv(16) Break	red/green	10	KBI-10304

Cat.# KBI-10305 RARA (17q21) Break

Acute Lymphoblastic Leukemia (ALL)

Acute lymphocytic leukemia, also called acute lymphoblastic leukemia, is a type of cancer that starts from white blood cells in the bone marrow. A number of recurring cytogenetic abnormalities are associated with distinct immunologic phenotypes of ALL and characteristic outcomes. The ETV6/AML1 fusion arising from the translocation t(12;21)(p13;q22) has been associated with a good outcome; the BCR/ ABL fusion of t(9;22)(g34;g11), rearrangements of the MLL gene (11q23), and abnormalities of the short arm of chromosomes 9 involving the tumor suppressor genes p16 (INK4A) have a poor prognosis.

ON TEL/AML t(12;21) Fusion

The t(12;21), a cryptic translocation rarely observed by conventional cytogenetics, was first identified by fluorescence in situ hybridization (FISH). In ALL blasts, this translocation fuses the 5' part of the TEL (ETV6) gene with almost the entire AML1 (CBFA2) gene, producing the chimeric transcript ETV6-CBFA2. The t(12;21) (p13;q22) has also been identified as the most frequent chromosomal abnormality in childhood ALL, affecting 20% to 25% of B-lineage cases.

The TEL/AML t(12;21) specific DNA probe is optimized to detect the reciprocal translocation t(12;21) (p13;q22) in a dual-color, dual-fusion assay.







TEL/AML t(12;21) Fusion probe hybridized to a normal metaphase (2R2G).

Literature[.]

Romana et al, 1995, Blood, 85: 3662-3670.

Ordering information	Color	Tests	Cat#
ON TEL/AML t(12;21) Fusion	red/green	10	KBI-10401

ON ETV6 (TEL) (12p13) Break

ETV6 (TEL) gene is the abbreviation for -ETS variant 6gene. It encodes an ETS family factor which functions as a transcriptional repressor in hematopoiesis and in vascular development. The gene is located on chromosome 12p13, and is frequently rearranged in human leukemias of myeloid or lymphoid origins. Also systematic deletion of the normal ETV6 allele in patients with ETV6-AML1 fusions can be found.

ON p16 (9p21) / 9q21

Hemizygous deletions and rearrangements of chromosome 9, band p21, are among the most frequent cytogenetic abnormalities detected in pediatric acute lymphoblastic leukemia (ALL). This deletion includes loss of the p16 (INK4A)/ p15 (INK4B) genes, which are cell cycle kinase inhibitors and important in leukemogenesis.

The p16 (9p21) specific DNA probe is optimized to detect copy numbers of the p16 (INK4A) gene region at region 9p21. The 9q21 region probe is included to facilitate chromosome identification.



ETV6 (TEL) (12p13) Break probe hybridized to patient material showing a translocation involving the ETV6 region at 12p13 (1RG1R1G). Image kindly provided by Magret Ratjen, Kiel

p16 (9p21) / 9q21 hybridized on patient material showing an isochromosome 9.

Image kindly provided by Dr Wenzel, Basel

Literature:

Golub et al, 1995, PNAS 92; 4917-4921 Ford et al, 2001, Blood 98; 558-564

Ordering information	Color	Tests	Cat#
ON ETV6 (TEL) (12p13) Break	red/green	10	KBI-10403

Literature:

Dreyling et al, 1995, Blood, 86: 1931-1938. Southgate et al, 1995, Br J Cancer, 72: 1214-1218.

Ordering information	Color	Tests	Cat#
ON p16 (9p21) / 9q21	red/green	10	KBI-10402

ON MLL (11q23) Break

Cat.# KBI-10303 MLL (11q23)

See description under AML on page 30.

Ordering information	Color	Tests	Cat#
ON MLL (11q23) Break	red/green	10	KBI-10303

Multiple Myeloma (MM)

The cytogenetic picture in multiple myeloma (MM) is highly complex, from which non-random numerical and structural chromosomal changes have been identified. Specifically, translocations involving the immunoglobulin heavy chain gene (IGH) at 14q32 and either monosomy or deletions of chromosome 13 have been observed in a significant number of patients. More recently several additional deletions or amplifications have been identified in MM which are currently investigated in large patient studies.

Note: Multiple Myeloma is a cancer of plasma cells. Analysis of such cells is hampered by their low frequency. Enrichment of plasma cells using CD138 is highly recommended.

ON BCR/ABL t(9;22)

The t(9;22) BCR/ABL translocation is present in about 5% of pediatric and up to 50% of adult ALL cases, and is associated with poor prognosis.

Cat.# KBI-10005 ON BCR/ABL t(9;22) Fusion Cat.# KBI-12005 ON BCR/ABL t(9;22) Fusion Cat.# KBI-10006 ON BCR/ABL t(9;22) TC, D-Fusion Cat.# KBI-10008 ON BCR/ABL t(9;22) DC, S-Fusion, ES Cat.# KBI-10009 ON BCR/ABL t(9;22) DC, S-Fusion See description under CML on page 13 and 14.

Ordering information	Color	Tests	Cat#
ON BCR/ABL t(9;22) Fusion	red/green	10	KBI-10005
ON BCR/ABL t(9;22) Fusion	red/green	20	KBI-12005
ON BCR/ABL t(9;22) TC, D-Fusion	red/green/blue	10	KBI-10006
ON BCR/ABL t(9;22) DC, S-Fusion, ES	red/green	10	KBI-10008
ON BCR/ABL t(9;22) DC, S-Fusion	red/green	10	KBI-10009

ON MYEOV/IGH t(11;14) Fusion

The most common chromosomal translocation in multiple myeloma (MM) is t(11;14), resulting in up-regulation of cyclin D1. In MM the breakpoints are scattered within a 360-kb region between CCND1 and MYEOV. This breakpoint is more proximal than the t(11;14) breakpoints observed in mantle cell lymphoma or other leukemias. Patients with MM who have t(11;14)(q13;q32) seem to have an aggressive clinical course.

The MYEOV/IGH t(11;14)(q13;q32) specific DNA probe is optimized to detect the reciprocal translocation t(11;14) in a dual-color, dual-fusion assay.

ON MAF/IGH t(14;16) Fusion

Chromosome translocations involving the immunoglobulin heavy chain gene (IgH) on 14q32 are a fundamental event in the pathogenesis of many B-cell malignancies. It often is preceded by a stable pre-malignant tumor called Monoclonal Gammopathy of Undetermined Significance (MGUS), which can sporadically progress to Multiple Myeloma (MM). One of the recurrent primary rearrangements involving the immunoglobulin heavy chain (IgH) locus on chromosome 14q32 identified in MGUS and MM tumors is the MAF/IgH t(14;16) translocation. Following MGUS appearance, the pathogenesis of multiple myeloma (MM) is thought to involve at least two pathways, which generate hyperdiploid (HRD) or nonhyperdiploid (NHRD) tumors, respectively.

The MAF/IGH is mainly present in NHRD tumors, providing important information on MM patient sub-types. Since these translocations are caused by aberrant IgH switch recombination, and possibly by aberrant somatic hypermutation in germinal center B cells, they provide information of an early and perhaps initiating event of transformation.





MYEOV/IGH t(11;14) Fusion probe hybridized to MM patient material showing t(11;14) translocation (2RG1R1G). Image kindly provided by Prof Jauch, Heidelberg

Cat.# KBI-10610 MAF/IGH t(14;16) Fusion





MAF/IGH t(14;16) Fusion probe hybridized to patient material showing a deletion of the MAF gene region at 16q23 (2R1G).

Literature:

Janssen et al., 2000, Blood, 95: 2691-2698. Fonseca et al, 2002, Blood, 99: 3735-3741

Ordering information	Color	Tests	Cat#
ON MYEOV/IGH t(11;14) Fusion	red/green	10	KBI-10605

Literature:

Chesi et at, 1998, Blood 91; 4457-4463 Sawyer et al, 1998, Blood 92; 4269-4278

Ordering information	Color	Tests	Cat#
ON MAF/IGH t(14;16) Fusion	red/green	10	KBI-10610

ON FGFR3/IGH t(4;14) Fusion

The t(4;14) translocation is undetectable by conventional cytogenetics. The breakpoints on chromosome 4 occur within an approximately 113-kb region located in small part of a conserved gene cluster including the transforming acidic coiled-coil protein 3 (TACC3), fibroblast growth factor receptor 3 (FGFR3), and multiple myeloma SET domain-containing protein (MMSET). The translocation predicts for poor survival and poor response to chemotherapy.

The FGFR3/IGH t(4;14)(p16;q32) specific DNA probe is optimized to detect the reciprocal translocation t(4;14) in a dual-color, dual-fusion assay

ON MM 11q23 / DLEU (13q14)

Hybridization results delineated 11q23 and 11q25 as the most frequently gained regions in MM, which supports a relevant pathogenetic role of genes in this region. Deletions of 13q14 are frequently detected by interphase FISH in patients with newly diagnosed MM, correlate with increased proliferative activity, and represent an independent adverse prognostic feature in MM.

The 11q23 specific DNA probe is optimized to detect copy numbers at 11q23. The DLEU (13q14) specific DNA region is optimized to detect copy numbers of the DLEU gene region at 13q14.

Cat.# KBI-10602 FGFR3/IGH t(4;14) Fusion



Cat.# KBI-10502 MM 11q23 / 13q14





FGFR3/IGH t(4;14) Fusion probe hybridized to MM patient material showing t(4;14) translocation (2RG1R1G). Image kindly provided by Prof Jauch, Heidelberg

Literature:

Chesi et al, 1997, Nat Genet, 16: 260-264. Finelli et al, 1999, Blood, 94: 724-732.

Ordering information	Color	Tests	Cat#
ON FGFR3/IGH t(4;14) Fusion	red/green	10	KBI-10602

Literature:

Gonzalez et al, 2004, Haematologica, 89: 1213-1218. Cremer et al, 2005, Genes Chrom Cancer, 44: 194-203.

Ordering information	Color	Tests	Cat#
ON MM 11q23 / DLEU (13q14)	red/green	10	KBI-10502

MM 11q23 / DLEU 13q14 probe hybridized to MM patient material

showing a 13q14 deletion (1R2G). Image kindly provided by Prof Jauch, Heidelberg

ON MM 1q21 / 8p21

Amplifications of 1q21 are concurrent with dysregulated expression of c-MAF, MMSET/FGFR3, or Del13 and represent an independent unfavorable prognostic factor. Allelic losses of the chromosome 8p21-22 have been reported as a frequent event in several cancers.

The 1q21 specific DNA probe is optimized to detect copy numbers at 1q21. The 8p21 specific DNA region is optimized to detect copy numbers at 8p21.

ON MM 19q13 / p53 (17p13)

P53 gene deletion, which can be identified by interphase FISH in almost a third of patients with newly diagnosed MM, is a novel prognostic factor predicting for short survival of MM patients treated with conventional-dose chemotherapy. Amplification of 19q13 has been reported in a variety of cancer.

The 19q13 specific DNA probe is optimized to detect copy numbers at 19q13. The p53 (17p13) specific DNA region is optimized to detect copy numbers of the p53 gene region at 17p13.





MM 1q21 / 8p21 hybridized to a normal metaphase (2R2G).

Literature:

Shaughnessy J., 2005, Hematology, 10 suppl 1: 117-126. Cremer et al, 2005, Genes Chrom Cancer, 44: 194-203.

Ordering information	Color	Tests	Cat#
ON MM 1q21 / 8p21	red/green	10	KBI-10503

Cat.# KBI-10509 MM 19q13 / p53 (17p13)





MM 19q13 / p53 (17p13) hybridized to patient material showing a p53 (17p13) deletion (1R2G). Image kindly provided by Prof Jauch, Heidelberg

Literature:

Drach et al, 1998, Blood, 92: 802-809. Cremer et al, 2005, Genes Chrom Cancer, 44: 194-203.

Ordering information	Color	Tests	Cat#
ON MM 19q13 / p53 (17p13)	red/green	10	KBI-10509

ON MM 15q22 / 6q21

Chromosome 6q amplifications encompassing 6q21-22 have been observed in MM including the same region as in CLL. Amplification including band 15q22 has been reported in MM.

The 15q22 specific DNA probe is optimized to detect copy numbers at 15q22. The 6q21 specific DNA region is optimized to detect copy numbers at 6q21.

ON MM 15q22 / 9q34

The hyperdiploid subtype in MM is defined by presence of multiple trisomic chromosomes. Combination of the chromosome 9q34 and 15q22 specific regions are important probe to recognize the hyperdiploid subtype in MM which is usually associated with a low frequency of IGH translocations.

The 15q22 and 9q34 probe is designed as a dual-color assay to detect amplifications at 15q22 and 9q34.

Cat.# KBI-10504 MM 15q22 / 6q21



Cat.# KBI-10508 MM 15q22 / 9q34





MM 15q22 / 6q21 hybridized to MM patient material with amplification of both critical regions 6q21 and 15q22. Image kindly provided by Prof Jauch, Heidelberg

Literature:

Cremer et al, 2005, Genes Chrom Cancer, 44: 194-203.

Ordering information	Color	Tests	Cat#
ON MM 15q22 / 6q21	red/green	10	KBI-10504



MM 15q22 / 9q34 hybridized to a normal interphase/ metaphase (2R2G)

Literature:

Cremer et al, 2005, Genes Chrom Cancer, 44: 194-203.

Ordering information	Color	Tests	Cat#
ON MM 15q22 / 9q34	red/green	10	KBI-10508

ON MM 1q21 / SRD (1p36)

Frequent loss of heterozygosity (LOH) on the short arm of chromosome 1 (1p) has been reported in a series of human malignancies. The combination with the potentially amplified 1g21 region allows to detect deletions at 1p36 and gain of 1q21 in a single FISH assay.

The 1q21 specific DNA probe is optimized to detect copy numbers at 1q21. The SRD 1p36 specific DNA Probe is optimized to detect copy numbers of 1p at region 1p36 containing the markers D1S2795 and D1S253.

ON IGH (14q32) Break

Multiple myeloma is characterized by complex rearrangements involving the IgH gene, particularly at the constant locus. The IgH rearrangement provides a useful marker of clonality in B-cell malignancies and amplification of this rearrangement is the method of choice to monitor the residual tumor cells in multiple myeloma.

The IGH (14q32) break probe is optimized to detect translocations involving the IGH gene region at 14g32 in a dual-color, split assay.



metaphase (2R2G).

Literature:

Cremer et al, 2005, Genes Chrom Cancer, 44: 194-203. Shaughnessy J., 2005, Hematology, 10 suppl 1: 117-126.

Ordering information	Color	Tests	Cat#
ON 1q21 / SRD (1p36)	red/green	10	KBI-10507

Literature:

Taniwaki et al, 1994, Blood, 83: 2962-1969. Moreau et al, 2002, Blood, 100: 1579-1583.

Ordering information	Color	Tests	Cat#
ON IGH (14q32) Break	red/green	10	KBI-10601

Lymphoma

Lymphoma is a type of cancer that originates in lymphocytes. Following WHO classification there are three large groups: B cell lymphomas, T cell and natural killer cell tumors, and Hodgkin lymphoma. The IGH gene locus at chromosome band 14q32 is very frequently involved in B-cell lymphoma.

ON MYC/IGH t(8;14) Fusion

The translocation t(8;14)(q24;q32) is the characteristic chromosomal aberration of Burkitt's-type of lymphomas. This translocation fuses the C-MYC gene at 8q24 next to the IGH locus at 14q32, resulting in overexpression of the transcription factor C-MYC. Detection of the t(8;14) is aimed to help in the diagnostic process of patients with highgrade B-cell lymphomas because treatment strategies differ between Burkitt and other high-grade lymphomas.

The MYC/IGH t(8;14)(q24;q32) specific DNA probe is optimized to detect the reciprocal translocation t(8;14) in a dual-color, dual-fusion assay.







MYC/IGH t(8;14) Fusion probe hybridized to a normal interphase/ metaphase (2R2G).

Literature:

Veronese et al, 1995, Blood, 85: 2132-2138. Siebert et al, 1998, Blood, 91: 984-990.

Ordering information	Color	Tests	Cat#
ON MYC/IGH t(8;14) Fusion	red/green	10	KBI-10603

ON BCL1/IGH t(11;14) Fusion

Mantle cell lymphoma is a subtype of non-Hodgkin lymphoma characterized by poor prognosis. Cytogenetically t(11;14) is associated with 75% of mantle cells lymphomas. The translocation breakpoints are scattered within the 120 kb BCL1 region adjacent to CCND1. The translocation leads to overexpression of cyclin D1 due to juxtaposition of the Ig heavy chain gene enhancer on 14q32 to the cyclin D1 gene on 11q13.

The BCL1/IGH t(11;14)(q13;q32) specific DNA probe is optimized to detect the reciprocal translocation t(11;14) in a dual-color, dual-fusion assay.

ON BCL2/IGH t(14;18) Fusion

The t(14;18) chromosomal translocation that results in the juxtaposition of the BCL2 proto-oncogene with the heavy chain JH locus. It a common cytogenetic abnormality in human lymphoma and is observed in about 85% of follicular lymphoma (FL) and up to one-third of diffuse lymphomas (DL). Two breakpoint regions cluster , major breakpoint region (mbr) within the 3' untranslated region of the BCL2 proto-oncogene accounts for approximately 60% of the cases and the minor cluster region (mcr) 30 kb 3' of BCL2 accounts for approximately 25% of the breakpoints.

The BCL2/IGH t(14;18)(q21;q32) specific DNA probe is optimized to detect the reciprocal translocation t(18;14) in a dual-color, dual-fusion assay.





Cat.# KBI-10606 BCL2/IGH t(14;18) Fusion





BCL1/IGH t(11;14) Fusion probe hybridized to a normal interphase/ metaphase (2R2G).

Literature:

Vaandrager et al, 1996, Blood, 88: 1177-1182.

Ordering information	Color	Tests	Cat#
ON BCL1/IGH t(11;14) Fusion	red/green	10	KBI-10604



BCL2/IGH t(14;18) probe hybridized to a normal interphase/ metaphase (2R2G).

Literature:

Poetsch et al, 1996, J Clin Oncol, 14: 963-969. Vaandrager et al, 2000, Genes Chrom Cancer, 27: 85-94.

Ordering information	Color	Tests	Cat#
ON BCL2/IGH t(14;18) Fusion	red/green	10	KBI-10606

ON IGH (14q32) Break

Chromosomal rearrangements involving the immunoglobulin heavy chain gene (IGH) at 14q32 are observed in 50% of patients with B-cell non-Hodgkin's lymphoma (NHL) and many other types of Lymphomas. More than 50 translocation partners with IGH have been described. In particular t(8;14), is associated with Burkitt's lymphoma, t(11;14) is associated with Mantle cell lymphoma, t(14;18) is observed in a high proportion of follicular lymphomas and t(3;14) is associated with Diffuse Large B-Cell Lymphoma.

The IGH (14q32) break probe is optimized to detect translocations involving the IGH gene region at 14q32 in a dual-color, split assay. Kreatech has developed this probe for the specific use on cell material (KBI-10601), or for the use on tissue (KBI-10729).

ON MALT (18q21) Break

Low grade malignant lymphomas arising from mucosa associated lymphoid tissue (MALT) represent a distinct clinicopathological entity. The three major translocations seen in MALT lymphomas are t(11;18)(q21;q21)/API2-MALT1, t(14;18)(q32;q21)/IGH-MALT1 and t(1;14)(p22;q32)/IGH-BCL10. A break or split probe for MALT (18q21) is best used to analyze translocation of the MALT gene on formalin fixed paraffin embedded tissue for routine clinical diagnosis.

Kreatech has developed this probe for the specific use on cell material (KBI-10608), or for the use on tissue (KBI-10731).





Cat.# KBI-10608 MALT (18q21) Break





MALT (18q21) Break probe hybridized to patient material showing a

translocation at 18q21 (1RG1RG).

IGH (14q32) Break probe hybridized to patient material showing a partial deletion of 14q32 (1RG1R).

Literature:

Taniwaki et al, 1994, Blood, 83: 2962-1969. Gozetti et al, 2002, Cancer Research, 62: 5523-5527.

Ordering information	Color	Tests	Cat#
ON IGH (14q32) Break	red/green	10	KBI-10601

Literature:

Morgan et al, 1999, Cancer Res, 59; 6205-6213 Dierlamm et al, 2000, Blood, 96; 2215-2218.

Ordering information	Color	Tests	Cat#
ON MALT (18q21) Break	red/green	10	KBI-10608

ON CCND1 (BCL1;11q13) Break

Besides the important functions in cellular growth, metabolism, and cellular differentiation, CCND1 (also known as Cyclin D1 or BCL1) can also function as a proto-oncogene, often dysregulated after re-arrangement by translocation. In fact, it can juxtapose into many different gene locus to drive tumorigenic effects. To date, the gene has been found to be rearranged in leukemias, in multiple myelomas (MM), and in some cases of benign parathyroid tumors. More specifically, the chromosomal translocation t(11;14)(q13:q32), involving rearrangement of the CCND1 locus, has been reported to be associated with human lymphoid neoplasia affecting mantle cell lymphomas (MCL). The rearrangement has been documented in 40-70% of cases of mantle cell lymphoma, whereas it only rarely occurs in other B cell lymphomas. In multiple myeloma, the same translocation t(11;14)(q13:q32) is the most common, with a reported frequency of 15% to 20% of the cases. For this reason, the CCND1 break apart probe KBI-10609 can be considered a very useful tool for routine diagnosis in MCL and Multiple myeloma, to be used in association to the related probes KBI-10604 and KBI-10605 probes that can detect more specifically the translocation t(11;14) in Mantle Cell Lymphoma (KBI-10604) and Multiple Myeloma (KBI-10605).



Cat.# KBI-10609 CCND1 (BCL1; 11q13) Break

CCND1 (BCL1; 11q13) Break probe hybridized to a normal metaphase (2R2G).

Literature:

Vaandrager et al, 1996, Blood, 88 (4); 1177-1182 Vaandrager et al, 1997, Blood, 89; 349-350

Ordering information	Color	Tests	Cat#
ON CCND1 (BCL1;11q13) Break	red/green	10	KBI-10609

ON BCL6 (3q27) Break

Chromosomal translocations involving band 3q27 with various different partner chromosomes represent a recurrent cytogenetic abnormality in B-cell non-Hodgkin's lymphoma. A FISH strategy using two differently labeled flanking BCL6 probes provides a robust, sensitive, and reproducible method for the detection of common and uncommon abnormalities of BCL6 gene in interphase nuclei. Kreatech has developed this probe for the specific use on cell material (KBI-10607), or for the use on tissue (KBI-10730).

Note: Two different breakpoint regions have been identified, the major breakpoint region (**MBR**) is located within the 5' noncoding region of the BCL6 proto-oncogene, while the atypical breakpoint region (**ABR**) is located approximately 250 kb proximal to the BCL6 gene. The Poseidon BCL6 (3q27) is designed in a way to flank both possible breakpoints, thereby providing clear split signals in either case.

Cat.# KBI-10607 BCL6 (3q27)



BCL6 (3q27) Break probe hybridized to patient material (1RG1R1G). Image kindly provided by Prof Siebert, Kiel

Literature:

Butler et al, 2002, Cancer Res, 62; 4089-4094. Sanchez-Izquierdo, 2001, Leukemia, 15; 1475-1484.

Ordering information	Color	Tests	Cat#
ON BCL6 (3q27) Break	red/green	10	KBI-10607

ON MYC (8q24) Break, TC

Rearrangements of the protooncogene C-myc (or MYC) have been consistently found in tumor cells of patients suffering from Burkitt's lymphoma. In cases with the common t (8;14) chromosomal translocation, the c-myc gene is translocated to chromosome 14 and rearranged with the immunoglobulin heavychain genes; the breakpoint occurs 5' to the c-myc gene and may disrupt the gene itself separating part of the first untranslated exon from the remaining two coding exons. In Burkitt's lymphoma showing the variant t (2;8) or t (8;22) translocations, the genes coding for the k and I mmunoglobulin light chain are translocated to chromosome 8. The rearrangement takes place 3' to the c-myc gene. At the present time the mechanism by which the oncogenic potential of the c-myc gene may be activated by these rearrangements is still controversial.

The MYC (8q24) break-apart probe is optimized to detect rearrangements involving the 8q24 locus in a triple-color, split assay on metaphase/interphase spreads, blood smears and bone marrow cells.



MYC (8q24) Break probe hybridized to patient material showing a 8q24 proximal break (1GBR1G1BR). Image kindly provided by Prof. Siebert, Kiel.

Literature:

Fabris et al, 2003, Genes Chromosomes Cancer 37 ; 261-269 Hummel et al., 2006, N Engl J Med 354 ;2419-30.

Ordering information	Color	Tests	Cat#
ON MYC (8q24) Break, TC	red/green	10	KBI-10611

ON FGFR3/IGH t(4;14) Fusion

Cat.# KBI-10602 FGFR3/IGH t(4;14) Fusion

See description under Multiple Myeloma on page 36.

Ordering information	Color	Tests	Cat#
ON FGFR3/IGH t(4;14) Fusion	red/green	10	KBI-10602

ON MYEOV/IGH t(11;14) Fusion

Cat.# KBI-10605 MYEOV/IGH t(11;14) Fusion

See description under Multiple Myeloma on page 35.

Ordering information	Color	Tests	Cat#
ON MYEOV/IGH t(11;14) Fusion	red/green	10	KBI-10605

ONCOLOGY PROBES - SOLID TUMORS





In solid tumors significantly high levels of chromosome abnormalities have been detected, but distinction between significant and irrelevant events has been a major challenge. Consequently, the application of cytogenetic technology as diagnostic, prognostic, or therapeutic tools for these malignancies has remained largely underappreciated. The emergence of FISH is particularly useful for solid malignancies, and the spectrum of their application is rapidly expanding to improve efficiency and sensitivity in cancer diagnosis, prognosis, and therapy selection, alone or with other diagnostic methods.

Description	Cat#	Page
ON ERBB2, Her-2/Neu (17q12) / SE 17	KBI-10701	47
ON EGFR, Her-1 (7p11) / SE 7	KBI-10702	59
ON CC hTERC (3q26) C-MYC (8q24) / SE 7 TC	KBI-10704	50
ON MYCN (2p24) / LAF (2q11)	KBI-10706	52
ON PPARγ (3p25) Break	KBI-10707	51
ON EWSR1 (22q12) Break	KBI-10708	55
ON hTERT (5p15) / 5q31 (tissue)	KBI-10709	59
ON p16 (9p21) / 9q21 (tissue)	KBI-10710	50
ON MLL (11q23) / SE 11	KBI-10711	53
ON SRD (1p36) / SE 1(1qh)	KBI-10712	53
ON SYT (18q11) Break	KBI-10713	56
ON CHOP (12q13) Break	KBI-10714	56
ON FUS (16p11) Break	KBI-10715	57
ON FKHR (13q14) Break	KBI-10716	57
ON MDM2 (12q15) / SE 12	KBI-10717	58
ON PTEN (10q23) / SE 10	KBI-10718	60

Breast Cancer

ON ERBB2, Her2/Neu (17q12) / SE 17

The Her2/Neu gene encodes a receptor tyrosine kinase involved in growth factor signaling. Overexpression of this gene is seen in about 20% of invasive breast cancers and is associated with poor survival. Her2 gene amplification is a permanent genetic change that results in this continuous overexpression of Her2. A drug called trastuzumab (commonly known as Herceptin) has been developed to be effective against Her2-positive breast cancer. Her2/Neu amplification is also observed in a variety of other tumors, such as prostate, lung, colon and ovary carcinoma.

The Her2/Neu (17q12) specific DNA probe is optimized to detect copy numbers of the Her2/Neu (ERBB2) gene region at region 17q12. The chromosome 17 satellite enumeration probe (SE 17) at D17Z1 is included to facilitate chromosome identification.

ON TOP2A (17q21) / SE 17

The Topoisomerase2A enzyme, which is vital for the cell because of its role in cell replication and repair, catalyzes the relaxation of supercoiled DNA molecules to create a reversible double-strand DNA break. This enzyme is also the target of a number of cytotoxic agents, namely TOP2A inhibitors (anthracyclines, etoposide, teniposide).

The dual-color probe is optimized to detect amplifications (copy numbers) or deletions of the TOP2A gene region at the 17q21. The chromosome 17 satellite enumeration probe (SE 17) at D17Z1 is included to facilitate chromosome identification.

Cat.# KBI-10701 ERBB2, Her2/Neu (17q12) / SE 17



ERBB2, Her2/Neu (17q12) / SE 17 probe hybridized to breast tumor tissue showing amplification of Her2/Neu (ERBB2)/ SE 17.

Literature:

Pauletti et al, 1996, Oncogene, 13: 63-72. Xing, et al, 1996, Breast Cancer Res Treat, 39: 203-212.

Ordering information	Color	Tests	Cat#
ON ERBB2, Her2/Neu (17q12) / SE 17	red/green	10	KBI-10701
ON ERBB2, Her2/Neu (17q12) / SE 17	red/green	50	KBI-14701

Cat.# KBI-10724 TOP2A (17q21) / SE 17



TOP2A (17q21) / SE 17 probe hybridized to breast tissue (2R2G).

Literature:

Järvinen et al, 1999, Genes, Chromosomes and Cancer 26; 142-150 Järvinen et al, 2000, Am. J. Pathology 156; 639-647

Ordering information	Color	Tests	Cat#
ON TOP2A (17q21) / SE 17	red/green	10	KBI-10724

ON TOP2A (17q21) / Her2 (17q12) / SE 17 Triple-Color Probe

The presence of both TOP2A amplification and deletion in advanced cancer are associated with decreased survival, and occur frequently and concurrently with Her2 gene amplification.

The TOP2A (17q21)/ Her2 (17q12)/ SE 17 probe is designed as a triple-color to detect amplification at 17q12 as well as amplifications or deletions at 17q21. The chromosome 17 satellite enumeration probe (SE 17) at D17Z1 in blue is included to facilitate chromosome identification.

ON ZNF217 (20q13) / 20q11

Zinc-finger protein 217 (ZNF217) is a Kruppel-like zinc-finger protein located at 20q13.2, within a region of recurrent maximal amplification in a variety of tumor types, and especially breast cancer cell lines and primary breast tumors. Copy number gains at 20q13 are also found in greater than 25% of cancers of the ovary, colon, head and neck, brain, and pancreas, often in association with aggressive tumor behavior. ZNF217 is considered a strong candidate oncogene that may have profound effects on cancer progression, which is transcribed in multiple normal tissues, and overexpressed in almost all cell lines and tumors in which it is amplified.

The ZNF217 (20q13) specific DNA probe is optimized to detect copy numbers of 20q at 20q13. The 20q11 probe is included to facilitate chromosome identification.

Cat.# KBI-10735 TOP2A (17q21) / Her2 (17q12) / SE 17



TOP2A (17q21)/ Her2(17q12) / SE 17 TC probe hybridized to breast tumor tissue showing amplification of TOP2A/Her2.

Literature:

Järvinen et al, 1999, Genes, Chromosomes and Cancer 26; 142-150 Järvinen et al, 2000, Am. J. Pathology 156; 639-647

Ordering information	Color	Tests	Cat#
ON TOP2A (17q21) / Her 2 / SE 17	red/green	10	KBI-10735

Cat.# KBI-10733 ZNF217 (20q13) / 20q11





ZNF217 (20q13) / 20q11 probe hybridized to tissue (2R2G).

Literature:

Tanner M et al, 2000, Clin Cancer Res, 6; 1833-1839 Ginestier C et al, 2006, Clin Cancer Res, 12; 4533-4544

Ordering information	Color	Tests	Cat#
ON ZNF217 (20q13) / 20q11	red/green	10	KBI-10733

ON FGFR1 (8p12) Break

FGFR1 (fibroblast growth factor receptor 1) expression has been shown to play pivotal roles in mammary development and breast cancer tumorigenesis. It has been shown that FGFR1 amplification is found in up to 10% of breast cancers and is significantly more prevalent in patients > 50 years of age and in tumors that lack HER2 expression. Even though the prognostic impact of FGFR1 amplification in breast cancer still remains unclear, the functional data demonstrating that FGFR1 signaling is required for the survival of breast cancer cells harboring FGFR1 amplification and the relatively high prevalence of FGFR1 amplification in breast cancer support the idea that this gene may be a useful therapeutic target for a subgroup of breast cancer patients with FGFR1 gene amplification. The FGFR1 (8p12) break-apart probe is optimized to detect translocations involving the FGFR1 gene region at 8p12 in a dual-color, split assay on metaphase/interphase spreads and paraffin embedded tissue sections.

Cat.# KBI-10737 FGFR1 (8p12) Break





FGFR1 (8p12) Break probe hybridized to patient material showing a break at 8p12 (1RG1R1G).

Literature:

Smedley et al, 1998, Hum Mol Genet. 7; 627-642. Sohal et al, 2001, Genes Chrom. Cancer 32; 155-163.

Ordering information	Color	Tests	Cat#
ON FGFR1 (8p12) Break	red/green	10	KBI-10737

Bladder Cancer

ON p16 (9p21) / 9q21 (tissue)

Homozygous and hemizygous deletions of 9p21 are the earliest and most common genetic alteration in bladder cancer. The p16 (INK4A) gene has been identified as tumor suppressor gene in this region which is commonly deleted in bladder cancer. The loss of DNA sequences on chromosomal bands 9p21-22 has been documented also in a variety of malignancies including leukemias, gliomas, lung cancers, and melanomas.

The p16 (9p21) specific DNA probe is optimized to detect copy numbers of the p16 gene region at region 9p21. The 9q21 region probe is included to facilitate chromosome identification.

Cat.# KBI-10710 p16 (9p21) / 9q21 (tissue)





p16 (9p21) / 9q21 (tissue) probe hybridized to tissue (2R2G).

Literature:

Stadler et al, 1994, Cancer Res, 54: 2060-2063. Williamson et al, 1995, Hum Mol Genet, 4: 1569-1577.

Ordering information	Color	Tests	Cat#
ON p16 (9p21) / 9q21 (tissue)	red/green	10	KBI-10710

Cervical Cancer

ON CC hTERC (3q26) / C-MYC (8q24) / SE 7 TC

Cervical cancer, a potentially preventable disease, remains the second most common malignancy in women worldwide. The most consistent chromosomal gain in aneuploid tumors of cervical squamous cell carcinoma mapped to chromosome arm 3q, including the human telomerase gene locus (hTERC) at 3q26. High-level copy number increases were also mapped to chromosome 8. Integration of HPV (Human Papilloma Virus) DNA sequences into C-MYC chromosomal regions have been repeatedly observed in cases of invasive genital carcinomas and in cervical cancers. The hTERC (3q26) specific DNA Probe is optimized to detect copy numbers of the hTERC gene region at region 3q26. The C-MYC (8q24) specific DNA probe is optimized to detect copy numbers of the C-MYC gene region at 8q24. The chromosome 7 satellite enumeration probe (SE 7) at D7Z1 is included as ploidy control.







CC hTERC (3q26) / C-MYC (8q24) / SE 7 TC probe hybridized to a PAP smear (destained) showing 3q26 and 8q24 amplification. The SE 7 control probe indicates a non-triploid karyotype (2B). Image kindly provided by Dr Weimer, Kiel

Literature:

Heselmeyer et al, 1996, PNAS, 93: 479-484. Herrick et al, 2005, Cancer Res, 65: 1174-1179.

Ordering information	Color	Tests	Cat#
ON CC hTERC (3q26) / C-MYC (8q24) / SE 7 TC	red/green	10	KBI-10704

Thyroid Carcinoma

Papillary thyroid carcinoma (PTC) is the most frequent primary carcinoma of the thyroid gland. The follicular carcinomas are associated with endemic goiter and a diet with low iodine intake. PTC, conversely, are multifocal and are associated with prior radiation and high iodine intake.

ON PPAR γ (3p25), Break

Follicular thyroid carcinoma is associated with the chromosomal translocation t(2;3)(q13;p25), fusing PAX8 (2q13) with the nuclear receptor, peroxisome proliferatoractivated receptor γ (PPAR γ). The close surrounding of PPAR is a breakpoint hot spot region, leading to recurrent alterations of this gene in thyroid tumors of follicular origin including carcinomas as well as adenomas with or without involvement of PAX8.

A break or split probe for PPAR γ is best used to analyze translocation of the PPAR γ (3p25) gene on formalin fixed paraffin embedded tissue for routine clinical diagnosis.

Cat.# KBI-10707 PPARy (3p25) Break





PPAR_γ (3p25) Break probe hybridized to patient material showing a translocation at 3p25 (1RG1R1G). Image kindly provided by Dr Valent, Paris

Literature:

French et al, 2003, Am J Pathol, 162; 1053-1060. Drieschner et l, 2006, Thyroid, 16; 1091-1096.

Ordering information	Color	Tests	Cat#
ON PPARγ (3p25), Break	red/green	10	KBI-10707

Neuroblastoma

According to the International Neuroblastoma Risk Grouping (INRG) Biology Committee MYCN remains the only genomic factor to be used currently for treatment stratification. Common data elements to be obtained by all groups include tumor cell ploidy and copy number/LOH status at chromosome bands 1p36, 11q23, and 17q23-25.

Literature:

Ambros et al, 2006, Advances in Neuroblastoma Research

ON MYCN (2p24) / LAF (2q11)

Amplification of the human N-myc protooncogene, MYCN, is frequently seen either in extrachromosomal double minutes or in homogeneously staining regions of aggressively growing neuroblastomas. MYCN amplification has been defined by the INRG as > 4-fold MYCN signals compared to 2q reference probe signals.

The MYCN (2p24) specific DNA probe is optimized to detect copy numbers of the MYCN gene region at 2p24. The LAF gene region probe at 2q11 is included to facilitate chromosome identification.

Cat.# KBI-10706 MYCN (2p24) / LAF (2q11)



MYCN (2p24) / LAF (2q11) hybridized to a cell line showing amplification of MYCN on chromosome 13 and 15. Image kindly provided by Pasteur Workshop 2008, Paris

Literature:

Shapiro et al, 1993, Am J Pathol, 142: 1339-1346. Corvi et al, 1994, PNAS, 91: 5523-5527..

Ordering information	Color	Tests	Cat#
ON MYCN (2p24) / LAF (2q11)	red/green	10	KBI-10706

ON SRD 1p36 / SE 1(1qh)

Neuroblastomas frequently have deletions of chromosome 1p and amplification of the N-myc oncogene. These deletions tend to be large and extend to the telomere, but a common region within sub-band 1p36.3 is consistently lost. Inactivation of a tumor suppressor gene within 1p36.3 is believed to be associated with an increased risk for disease relapse. The 1p36 specific DNA probe has recently been changed to cover the recently described smallest region of consistent deletion **(SRD)** between D1S2795 and D1S253.

The SRD (1p36) specific DNA probe is optimized to detect copy numbers of the 1p36 region on chromosome 1. The chromosome 1 satellite enumeration probe (SE 1) at 1qh is included to facilitate chromosome identification.

ON MLL (11q23) / SE 11

Deletions of the long arm of chromosome 11 (11q) have been noted in primary neuroblastomas. It is assumed that a tumor suppressor gene mapping within 11q23.3 is commonly inactivated during the malignant evolution of a large subset of neuroblastomas, especially those with unamplified MYCN.

The MLL (11q23) specific DNA probe is optimized to detect amplification or deletion involving the MLL gene region at 11q23 in a dual-color assay on metaphase/interphase spreads, blood smears and bone marrow cells. The Chromosome 11 Satellite Enumeration probe (SE 11) at D11Z1 is included to facilitate chromosome identification.





SRD (1p36) / SE 1 probe hybridized to a normal metaphase (2R2G).

Literature:

Caron et al, 1993, Nat Genet, 4: 187-190. Cheng et al, 1995, Oncogene, 10: 291-297. White et al, 2005, Oncogene, 24: 2684-2694

Ordering information	Color	Tests	Cat#
ON SRD (1p36) / SE1(1qh)	red/green	10	KBI-10712

Cat.# KBI-10711 MLL (11q23) / SE 11





MLL (11q23) / SE 11 hybridized to normal interphases (2R2G).

Literature:

Guo et al, 1999, Oncogene, 18: 4948-4957. Maris et al, 2001, Med Pediatr Oncol, 36: 24-27.

Ordering information	Color	Tests	Cat#
ON MLL (11q23) / SE11	red/green	10	KBI-10711

ON p53 (17p13) / MPO (17q22) "ISO 17q"

Gain of genetic material from chromosome arm 17q (gain of segment 17q21-qter) is the most frequent cytogenetic abnormality of neuroblastoma cells. In multivariate analysis, 17q gain was more strongly associated with adverse outcome than was either stage (Stage 4 vs other combined) or 1p status.

ON MDM4 (1q32) / SE1

MDM4 (also known as MDMX, murine double minute gene) is a relative of MDM2 that was identified on the basis of its ability to physically interact with p53. MDM4, like MDM2, acts as a key negative regulator of p53 by interfering with its transcriptional activity. MDM4 amplification and/or overexpression occurs in several diverse tumors. Studies showed an increased MDM4 copy number in 65% of human retinoblastomas compared to other tumors, qualifying MDM4 as a specific chemotherapeutic target for treatment of this tumor.

The MDM4 (1q32) specific DNA probe is designed as a dualcolor assay to detect amplification at 1q32. The chromosome 1 Satellite Enumeration (SE 1) probe at 1qh is included to facilitate chromosome identification.

Cat.# KBI-10011 p53 (17p13) / MPO (17q23) "ISO 17q"

See probe description under CML on page 15.

Cat.# KBI-10736 MDM4 (1q32) / SE 1



MDM4 (1q32) / SE 1 probe hybridized to paraffin embedded tissue (2R2G).

Literature:

Bown et al, 1999, N Engl J Med, 340: 1954-1961. O'Neill, 2001, Genes Chrom Cancer, 30: 87-90.

Ordering information	Color	Tests	Cat#
ON p53 (17p13) / MPO (17q22) "ISO 17q"	red/green	10	KBI-10011

Literature:

Riemenschneider et al, 1999, Cancer Res. 59 ; 6091-6096 Danovi et al, 2004, Mol.Cell.Bio. 24; 5835-5843

Ordering information	Color	Tests	Cat#
ON MDM4 (1q32) / SE 1	red/green	10	KBI-10736

Sarcoma

Sarcoma is a general class of less common cancers in which the cancer cells arise from or resemble normal cells in the body known as "connective tissues" (fat, muscle, blood vessels, deep skin tissues, nerves, bones, and cartilage). The benign and malignant forms have related karyotypic changes which provide an important resource for identifying the additional genetic changes that occur in the malignant compared with the benign form. In fact, the molecular biology of soft-tissue sarcomas has provided the perfect example of how cytogenetic and molecular approaches can contribute toward a clearer understanding of the development of soft-tissue sarcomas.

ON EWSR1 (22q12) Break

The translocation t(11;22)(q24;q12) has been independently described by several groups in the malignant cells of patients with Ewing sarcoma and can be detected in more than 90% of these tumors. The translocation involves the fusion of the human FLI1 gene on chromosome 11, with the coding sequence of the EWS gene on chromosome 22. In about 5% of the cases, however, the EWS gene is involved in variant translocations. These translocations, t(21;22)(q12;q12) and t(7;22)(p22;q12) result in the fusion of EWS with ERG and ETV1, respectively.

The Poseidon break or split probe for EWS is optimized to identify translocation involving the EWS gene on formalin fixed paraffin embedded tissue for routine clinical diagnosis.

Cat.# KBI-10708 EWS (22q12) Break





EWS (22q12) Break probe hybridized to patient material showing a translocation involving the EWRS1 gene region at 22q12 (1RG1R1G).

Literature:

Zucman-Rossi, et al, 1998, PNAS, 95; 11786-11791. Bernstein et al, 2006, Oncologist, 11; 503-519.

Ordering information	Color	Tests	Cat#
ON EWSR1 (22q12) Break	red/green	10	KBI-10708

ON SYT (18q11) Break

The characteristic chromosomal abnormality in synovial sarcoma is t(X;18)(p11.2;q11.2) present in 90% of patients. This translocation results in the fusion of the chromosome 18 SYT gene to either of two distinct genes, SSX1 or SSX2, located on the X chromosome.

A break or split probe for SYT is best used to analyze translocation of the SYT (SS18) gene on formalin fixed paraffin embedded tissue for routine clinical diagnosis.

ON CHOP (12q13) Break

Liposarcoma is one of the most frequent sarcomas in adults, representing 10 to 16 percent of soft tissue sarcomas. Most patients with round cell/myxoid liposarcoma have an acquired t(12;16)(CHOP-FUS) or t(12;22)(CHOP-EWS) translocation, both of which involve the CHOP gene at 12q13. A break or split probe for CHOP is best used to analyze translocation of the CHOP (12q13) gene on formalin fixed paraffin embedded tissue for routine clinical diagnosis.

The CHOP (12q13) Break probe is optimized to detect translocations involving the CHOP gene region at 12q13 in a dual-color, split assay.



Cat.# KBI-10714 CHOP (12q13) Break





CHOP (12q13) Break probe hybridized to a normal metaphase (2RG).

Literature:

Kawai et al, 1998, NEJM, 338; 153-160 Surace et al, 2004, LabInvest., 84; 1185-1192

Ordering information	Color	Tests	Cat#
ON SYT (18q11) Break	red/green	10	KBI-10713

Literature:

Panagopoulos et a, 1994, Cancer Res, 54; 6500-6503. Schoenmakers et al, 1994, Genomics, 20; 210-222.

Ordering information	Color	Tests	Cat#
ON CHOP (12q13) Break	red/green	10	KBI-10714

ON FUS (16p11) Break

The FUS gene was originally shown to be rearranged in myxoid liposarcomas harboring a t(12;16)(q13;p11) translocation. FUS has been shown to be involved also in other recombinations: with ERG in acute myeloid leukemia carrying a t(16;21), with ATF1 in band 12q13 in angiomatoid fibrous histiocytoma, and with CREB3L2 in fibromyxoid sarcoma.

A break or split probe for FUS is best used to analyze translocation of the FUS (16p11) gene on formalin fixed paraffin embedded tissue for routine clinical diagnosis.

ON FKHR (13q14) Break

The t(2;13) is associated with alveolar rhabdomyosarcomas. This translocation results in the formation of a chimeric transcript consisting of the 5' portion of PAX3, including an intact DNA-binding domain fused to the FKHR gene on chromosome 13. The t(1;13)(p36;q14) also seen in alveolar rhabdomyosarcomas results in the fusion of another member of the PAX family, PAX7 to the FKHR gene on chromosome 13.

A break or split probe for FKHR is best used to analyze translocation of the FKHR (13q14) gene on formalin fixed paraffin embedded tissue for routine clinical diagnosis.



Shing et al, 2003, Cancer Res, 63: 4568-4576. Storlazzi et al, 2003, Hum. Mol. Genet., 12: 2349-2358.

Ordering information	Color	Tests	Cat#
ON FUS (16p11), Break	red/green	10	KBI-10715

Barr et al, 1996, Hum, Mol. Genet., 5: 15-21. Coignet et al, 1999, Genes Chrom. Cancer, 25; 222-229.

Ordering information	Color	Tests	Cat#
ON FKHR (13q14) Break	red/green	10	KBI-10716

Incology Probes - Solid Tumors

57

ON MDM2 (12q15) / SE 12

Fibrosarcoma is a rare soft-tissue tumor composed of fascicles of spindled fibroblast-like cells. Gains and high-level amplifications of 12q14–22 were the most common genomic imbalances, and reflected MDM2 amplification, thereby indicating the importance of this gene in the evolution of fibrosarcomas.

The MDM2 (12q15) specific DNA probe is optimized to detect copy numbers of the MDM2 gene region at region 12q15. The Chromosome 12 Satellite Enumeration probe (SE 12) at D12Z3 is included to facilitate chromosome identification.

ON CDK4 (12q13) / SE 12

Amplification of the CDK4 gene region at 12q13-q15 has been observed in several types of cancer, especially in gliomas and sarcomas. CDK4 codes for a cyclin dependent kinase which is involved in controlling progression through the G1 phase of the cell cycle. The oncogenic potential of CDK4 activation has been related to the deregulation of the G1 phase by increasing the hyperphosphorylation of retinoblastoma tumor suppressor protein helping to cancel its growth-inhibitory effects.

The CDK4 (12q13) specific DNA probe is optimized to detect copy numbers of the CDK4 gene region at 12q13. The chromosome 12 satellite enumeration probe (SE 12) at D12Z3 is included to facilitate chromosome identification.

Cat.# KBI-10717 MDM2 (12q15) / SE 12



Cat.# KBI-10725 CDK4 (12q13) / SE 12





MDM2 (12q15) / SE 12 hybridized to a normal metaphase (2R2G).

Literature:

Mitchell et al, 1995, Chrom. Res., 3; 261-262. Reifenberger et al, 1996, Cancer Res., 15; 5141-5145.

Ordering information	Color	Tests	Cat#
ON MDM2 (12q15) / SE12	red/green	10	KBI-10717



CDK4 (12q13) / SE 12 probe hybridized to liposarcoma tissue showing multiple amplification involving the CDK4 gene region at 12q13 (3+R2G). Image kindly provided by Dr. Sapi, Hungary.

Literature:

Kuhnen et al, 2002, Virchows Arch 441 ; 299-302 Shimada et al, 2006, Hum Path 37(9) ; 1123-1129

Ordering information	Color	Tests	Cat#
ON CDK4 (12q13) / SE 12	red/green	10	KBI-10725

Different Cancer types

ON EGFR, Her-1 (7p11) / SE 7

Epidermal growth factor receptor (EGFR) is a cell membrane protein, providing signal transduction and cell growth. It is a member of the Her or Erb-B family of type I receptor tyrosine kinases and implicated in the development and progression of cancer of the lung, breast, intestine, and other organs. EGFR was found to act as a strong prognostic indicator in head and neck, ovarian, cervical, bladder and oesophageal cancers. In these cancers, increased EGFR expression was associated with reduced recurrence-free or overall survival.

The EGFR (7p11) specific DNA probe is optimized to detect copy numbers of the EGFR (Her1) gene region at region 7p11. The chromosome 7 satellite enumeration probe (SE 7) at D7Z1 is included to facilitate chromosome identification.

Cat.# KBI-10702 EGFR, Her1 (7p11) / SE 7





EGFR, Her1 (7p11) / SE 7 hybridized to colon tissue (2R2G).

Literature:

Wang et al, 1993, Jpn J Hum Genet, 38: 399-406. Nicholoson et al, 2001, Eur J Cancer, 37: 9-15.

Ordering information	Color	Tests	Cat#
ON EGFR, Her-1 (7p11) / SE 7	red/green	10	KBI-10702

ON hTERT (5p15) / 5q31 (tissue)

Amplification of the hTERT gene at band 5p15 has been observed in a variety of cancer, particularly lung cancer, cervical tumors, and breast carcinomas. Several studies have revealed a high frequency of hTERT gene amplification in human tumors, which indicates that the hTERT gene may be a target for amplification during the transformation of human malignancies and that this genetic event probably contributes to a dysregulation of hTERT/telomerase occurring in a subset of human tumors.

The hTERT (5p15) probe is designed as a dual-color assay to detect amplification at 5p15. The CDC25C/EGR1 (5q31) gene region probe is included as internal control.

Cat.# KBI-10709 hTERT (5p15) / 5q31 (tissue)



hTERT (5p15) / 5q31 (tissue) probe hybridized to paraffine embedded tissue (2R2G).

Literature:

Bryce et al, 2000, Neoplasia, 2; 197-201. Zhang et al, 2000, Cancer Res, 60; 6230-6235.

Ordering information	Color	Tests	Cat#
ON hTERT (5p15) / 5q31 (tissue)	red/green	10	KBI-10709

ON PTEN (10q23) / SE 10

The gene 'phosphatase and tensin homolog deleted on chromosome 10' (PTEN), is a tumor suppressor located at chromosome 10q23, that plays an essential role in the maintenance of chromosomal stability, cell survival and proliferation. Loss of PTEN has been found in a wide number of tumors, and his important role is demonstrated by the fact that it is the second most frequently mutated gene after p53. Loss of PTEN significantly correlates with the advanced forms of gliomas, but also of prostate cancer and breast tumors.

The PTEN (10q23) specific DNA probe is optimized to detect copy numbers of the PTEN gene region at region 10q23. The Chromosome 10 Satellite enumeration probe (SE 10) at D10Z1 is included to facilitate chromosome identification.

ON C-MET (7q31) / SE 7

The C-MET proto-oncogene is a receptor-like tyrosine kinase that drives a physiological cellular program important for development, cell movement, cell repair, cellular growth. Aberrant execution of the program has been associated to neoplastic transformation, invasion and metastasis. Activation of C-MET has been reported in a significant percentage of human cancers and is amplified during the transition between primary tumors and metastasis.

The C-MET (7q31) specific DNA probe is optimized to detect copy numbers of the C-MET gene region at region 7q31. The Chromosome 7 Satellite enumeration probe (SE 7) at D7Z1 is included to facilitate chromosome identification.







PTEN (10q23) / SE 10 probe hybridized to prostate cancer material showing deletion of PTEN gene region at 10q23 (1R2G). Image kindly provided by Portuguese Cancer Inst., Porto

Literature:

Cairns et al, 1997, Cancer Res, 57 ; 4997-5000 Hermans et al, 2004, Genes Chrom Cancer, 39; 171-184

Ordering information	Color	Tests	Cat#
ON PTEN (10q23) / SE 10	red/green	10	KBI-10718

Cat.# KBI-10719 C-MET (7q31) / SE 7



C-MET (7q31) / SE 7 probe hybridized to colon tissue (2R2G).

Literature:

Hara et al, 1998, Lab Invest 78; 1143-1153. Tsugawa et al., 1998, Oncology 55; 475-481.

Ordering information	Color	Tests	Cat#
ON C-MET (7q31) / SE 7	red/green	10	KBI-10719

ON AURKA (20q13) / 20q11

Aurora kinase A (AURKA) covers the fundamental role of regulating proper centrosome function, important to maintain genomic stability during cell division and to ensure equal segregation of replicated chromosomes to daughter cells. Deregulated duplication and distribution of centrosomes has been implicated in mechanisms leading to mitotic spindle aberrations, aneuploidy, and genomic instability that are seen in many different tumor types. Consistent with this, AURKA amplification has been detected in approximately 12% of primary breast tumors, as well as in breast, ovarian, colon, prostate, neuroblastoma and cervical cancer cell lines. Recent investigations on new drugs developments have focused on the importance of aurora kinases for tumor suppression.

The AURKA (20q13) specific DNA probe is optimized to detect copy numbers of the AURKA gene region at region 20q13. The 20q11 specific DNA probe is included to facilitate chromosome identification.

ON AURKB (17p13) / SE 17

Aurora kinase B (AURKB) localizes to microtubules, and is a key regulator of the mitotic cell division and chromosome segregation processes. Gain of function of AURKB correlates with cell proliferation, induction of multinuclear cells, and chromosomal instability. The significant interest of the gene in cancer diagnostics is related to the driving function of AURKB in tumor progression, histological differentiation, and metastasis. AURKB is predictive for the aggressive recurrence of many different types of tumors, including hepatocellular carcinoma and oral squamous cell carcinoma. Recently new drugs have been under investigation for their capacity of interfering with the aurora kinases activity related to tumor-suppressor effects.

The AURKB (17p13) specific DNA probe is optimized to detect copy numbers of the AURKB gene region at region 17p13. The Chromosome 17 Satellite Enumeration (SE 17) probe at D17Z1 is included to facilitate chromosome identification.

Color

red/areen

Tests Cat#

10

KBI-10722

Ordering information

ON AURKB (17p13) / SE 17



Ordering information	Color	Tests	Cat#
ON AURKA (20q13) / 20q11	red/green	10	KBI-10721

ON CCND1 (11q13) / SE 11

CCND1 (also named Cyclin D1 or BCL1) is a key cell cycle regulator of the G1 to S phase progression. The binding of cyclin D1 to cyclin-dependent kinase (CDKs) leads to the phosphorylation of retinoblastoma protein (pRb), subsequently triggering the release of E2F transcription factors to allow G1 to S phase progression of the cell cycle. Consistent with this function, overexpression of cyclin D1 results in a more rapid progression from the G1 to S phase transition and in a reduced serum dependency in fibroblast cells, characteristics typically seen in cancer cells. Amplification of cyclin D1 plays pivotal roles in the development of a subset of human cancers including parathyroid adenoma, breast cancer, colon cancer, lymphoma, melanoma, and prostate cancer.

The CCND1 (11q13) specific DNA Probe is optimized to detect copy numbers of the CCND1 gene region at region 11q13. The Chromosome 11 Satellite Enumeration (SE 11) probe at D11Z1 is included to facilitate chromosome identification.

Cat.# KBI-10734 CCND1 (11q13) / SE 11





CCDN1 (11q13) / SE 11 probe hybridized to patient interphases/ metaphase showing CCDN1 (11q13) amplification with polyploidy for chromosome 11.

Literature:

Okami et al, 1999, Oncogene 18; 3541-3545. Freier et al, 2003, Cancer Res; 1179-1182

Ordering information	Color	Tests	Cat#
ON CCND1 (11q13) / SE 11	red/green	10	KBI-10734

TMPRSS2-ERG (21q22) Del, Break, TC

The transmembrane protease serine 2 gene (TMPRSS2) is involved in gene fusions with ERG, ETV1 or ETV4 in prostate cancer. In recent studies it has been described that the expression of the TMPRSS2-ERG fusion gene is a strong prognostic factor for the risk of prostate cancer recurrence in prostate cancer patients treated by surgery.

The TMPRSS2-ERG rearrangement probe is optimized to detect the deletion between TMPRSS2 and ERG at 21q22 associated with the TMPRSS2-ERG fusion in a triple-color deletion assay. It also detects translocations involving the TMPRSS2 region such as ETV1 t(7;21), or ETV4 t(17;21).

ON IGH (14q32) Break (tissue)

Chromosomal rearrangements involving the immunoglobulin heavy chain gene (IGH) at 14q32 are observed in 50% of patients with B-cell non-Hodgkin's lymphoma (NHL) and many other types of Lymphomas. More than 50 translocation partners with IGH have been described. In particular t(8;14), is associated with Burkitt's lymphoma, t(11;14) is associated with Mantle cell lymphoma, t(14;18) is observed in a high proportion of follicular lymphomas and t(3;14) is associated with Diffuse Large B-Cell Lymphoma.

The IGH (14q32) break probe is optimized to detect translocations involving the IGH gene region at 14q32 in a dual-color, split assay. Kreatech has developed this probe for the specific use on cell material (KBI-10601), or for the use on tissue (KBI-10729).

Cat.# KBI-10726 TMPRSS2-ERG (21q22) Del, Break, TC



TMPRSS2-ERG (21q22) rearrangement probe hybridized to prostate carcinoma tissue showing a deletion of the TMPRSS2 (21q22) gene region associated with TMPRSS2-ERG fusion (1RGB 1RB). Image kindly provided by Dr. Texeira, Porto.

Literature:

Perner et al, 2006 Cancer Res 66; 8337-8341 Hermans et al, 2006, Cancer Res 66; 10658-10663 Attard et al, 2008, Oncogene 27; 253-263

Ordering information	Color	Tests	Cat#
TMPRSS2-ERG (21q22) Del, Break, TC	red/green	10	KBI-10726

Cat.# KBI-10729 ON IGH (14q32) Break (tissue)





IGH (14q32) Break probe hybridized to patient material showing a partial deletion of 14q32 (1RG1R).

Literature:

Taniwaki et al, 1994, Blood, 83: 2962-1969. Gozetti et al, 2002, Cancer Research, 62: 5523-5527.

Ordering information	Color	Tests	Cat#
ON IGH (14q32) Break (tissue)	red/green	10	KBI-10729

ON BCL6 (3q27) Break (tissue)

Chromosomal translocations involving band 3q27 with various different partner chromosomes represent a recurrent cytogenetic abnormality in B-cell non-Hodgkin's lymphoma. A FISH strategy using two differently labeled flanking BCL6 probes provides a robust, sensitive, and reproducible method for the detection of common and uncommon abnormalities of BCL6 gene in interphase nuclei. Kreatech has developed this probe for the specific use on cell material (KBI-10607), or for the use on tissue (KBI-10730).

Note: Two different breakpoint regions have been identified, the major breakpoint region (**MBR**) is located within the 5' noncoding region of the BCL6 proto-oncogene, while the atypical breakpoint region (**ABR**) is located approximately 250 kb proximal to the BCL6 gene. The Poseidon BCL6 (3q27) is designed in a way to flank both possible breakpoints, thereby providing clear split signals in either case.

ON MALT (18q21) Break (tissue)

Low grade malignant lymphomas arising from mucosa associated lymphoid tissue (MALT) represent a distinct clinicopathological entity. The three major translocations seen in MALT lymphomas are t(11;18)(q21;q21)/API2-MALT1, t(14;18)(q32;q21)/IGH-MALT1 and t(1;14)(p22;q32)/ IGH-BCL10. A break or split probe for MALT (18q21) is best used to analyze translocation of the MALT gene on formalin fixed paraffin embedded tissue for routine clinical diagnosis.

Kreatech has developed this probe for the specific use on cell material (KBI-10608), or for the use on tissue (KBI-10731).

Cat.# KBI-10730 ON BCL6 (3q27) Break (tissue)





BCL6 (3q27) Break probe hybridized to patient material showing both normal (2RG) and aberrant signals (1RG1R1G). Image kindly provided by Prof Siebert, Kiel

Literature:

Butler et al, 2002, Cancer Res, 62; 4089-4094. Sanchez-Izquierdo, 2001, Leukemia, 15; 1475-1484.

Ordering information	Color	Tests	Cat#
ON BCL6 (3q27) Break (tissue)	red/green	10	KBI-10730

Cat.# KBI-10731 ON MALT (18q21) Break (tissue)





MALT (18q21) Break tissue probe hybridized to paraffin embedded material (2RG).

Literature:

Morgan et al, 1999, Cancer Res, 59; 6205-6213 Dierlamm et al, 2000, Blood, 96; 2215-2218.

Ordering information	Color	Tests	Cat#
ON MALT (18q21) Break (tissue)	red/green	10	KBI-10731

ON p53 (17p13) / SE 17 (tissue)

The p53 tumor suppressor gene at 17p13, has been shown to be implicated in the control of normal cellular proliferation, differentiation, and apoptosis. Allelic loss, usually by deletion, and inactivation of p53 have been reported in numerous tumor types and are associated with poor prognosis in CLL.

The p53 (17p13) specific DNA probe is optimized to detect copy numbers of the p53 gene region at 17p13. The chromosome 17 satellite enumeration probe (SE 17) at D17Z1 is included to facilitate chromosome identification.

ON TFE3 (Xp11) Break

Abnormalities of Xp11.2 region have often been observed in papillary renal cell carcinomas and are sometimes the sole cytogenetic abnormality present. The transcription factor binding to IGHM enhancer 3 (TFE3) gene, which encodes a member of the helix-loop-helix family of transcription factors, is located in this critical region and can be fused to various other chromosomal regions by translocation. Known fusion partners are NONO (Xq12), PRCC (1q21), SFPQ (1p34), CLTC (17q23) and ASPSCR1 (17q25).

The TFE3 (Xp11) Break probe is optimized to detect translocations involving the TFE3 gene region at Xp11.2 in a dual-color, split assay.

